- GB** Timesco Single Use Blade & Handle Instructions for Use

A laryngoscope is used to examine and visualise a patient's upper airway and aid placement of a tracheal tube. Only to be used by a qualified physician. (N.B. USA and Canadian Federal Law restricts this device to sale by or on the order of a physician.)

These blades and handles conform to the ISO 7376:2020 fitting, only use with systems that meet this standard. Ensure you use compatible systems.

- Supplied non-sterile, packed in clean conditions.
- Do not reuse, sterilise or reprocess.

Pre-use checks – before use on a patient

- Visually inspect the product for any sharp edges, loose parts or damage.
- Fit the blade to the handle and check that the light functions correctly.

Cautions

- This device is for single product and patient use only. DO NOT RE-USE OR REPROCESS. Re-use of a single use device carries with it the potential risk of contamination of the device and/or risk of patient/ user infection or cross infection including but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination, cross-contamination and/or cross-infection may lead to injury, illness, or death of the patient. Reprocessing of the device may not be effective and may compromise the structural integrity of the device and/or lead to device failure, which in turn may lead to patient injury, illness or death. The manufacturer shall not be liable for any damages caused by re-use or reprocessing of this device or accessories.
- Bulbs left illuminated in an exposed position may generate heat sufficient to burn human tissue.
- Do not use in the vicinity of MRI scanners or intense magnetic fields.
- Repeated testing before use may reduce the life of the product or result in operational failure.
- Incorrectly fitting the blade to the handle may result in a non-functioning or damaged blade and/or handle.
- Always keep spare blades and handles available in case of failure.

Fitting & removing a blade from a handle


- Hook the blade onto the hook bar on the handle (1)
- Push the back of the blade block until the blade 'clicks' securely onto the hook bar (2)
- Lift the tip of the blade until it locks into the operating position and the bulb illuminates (3).

To remove the blade, reverse the instructions above.

Fitting batteries in handles that DO NOT CONTAIN PRELOADED BATTERIES: Use only Alkaline batteries. Use the correct size and quantity of batteries per handle.

Used blades, handles and batteries should be discarded as clinical waste in accordance with local policies.

Please report any serious incidents caused by this product to Timesco and the national regulatory body.


	See instructions for use	MD	Medical device
	Non-Sterile	REF	Catalogue number
	Do not re-use	LOT	Batch code
	Use by date (where applicable)	UDI	Unique device identifier
	Latex free		Type BF applied part
	Manufacturing date (where applicable)		Keep dry
	Manufacturer		Not to be disposed of with household waste

- Do not use if packaging has been opened or damaged

Rx ONLY Federal (U.S.A) law restricts this device to sale by or on the order of a physician

EC REP Authorized Representative in the European Community

	Importer
---	----------

Recommended operating environment		
Operating temperature and humidity:		
	10°C to 40°C	30% to 70%
Storage/Transport temperature and humidity:		
	-20°C to 60°C	10% to 70%
Atmospheric pressure:		
	70kPa to 106kPa	

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