





# GENERAL MEDICAL

PROFESSIONAL ESSENTIALS  
TIMESCO

## PRODUCT RANGE

## PAGE

Sphygmomanometers	2-4	
Stethoscopes	5-8	
Sphygmomanometer & Stethoscopes Sets & Pink Products	9	
Pen Torches	10	
Watches	11	
Scissors	12	
Tourniquets	13	
Suction Devices	14	
Resuscitation Aids	14-15	

## WHO WE ARE



### REPUTABLE BRAND.

As a family owned business, Timesco have been at the forefront of medical product design & innovation since 1964 & are recognised as one of the most progressive & innovative companies within the industry. We are proud to manufacture & supply an unrivalled portfolio of medical products to medical professionals in Anaesthesia, Surgery, Podiatry & General Medical.



### END USER SATISFACTION.

It's important that the end user is satisfied with the products that we supply. From concept through to design & manufacture, our aim is to deliver products made to exceed internationally recognised quality standards & customer expectations.



### EFFICIENT LOGISTICS.

All our products are available for direct orders, or we can arrange delivery through our network of global distributors. Many of our General Medical products are available for procurement through various framework agreements. You can find the corresponding product codes in this brochure.



### SALES SUPPORT.

We have a dedicated external field sales team strategically located, offering extensive global coverage (Territory Manager details are available on request).

We also have an internal sales customer support team based at our head office.

Our sales and support team are available Monday to Friday from 9am GMT to 5pm GMT.

Tel: 01268 297710 | Email: [sales@timesco.com](mailto:sales@timesco.com) | Live chat: [www.timesco.com](http://www.timesco.com)



## EMERALD CLIP-ON ANEROID SPHYGMOMANOMETER

- 2 tube system
- Separate bellows & gauge tubes
- Spring loaded clip-on rear of dial for one-handed use
- 48mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult black nylon calibrated cuff
- Zip-up pouch for convenient storage
- 1 year guarantee

DESCRIPTION	PRODUCT CODE
Emerald Sphygmomanometer	D05.105



## SAPPHIRE ANEROID SPHYGMOMANOMETER

- 1 tube system
- Combined bellows & gauge tubes
- Easy one-handed use
- 48mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult nylon calibrated cuff
- Zip-up pouch for convenient storage
- 1 year guarantee

COLOUR	PRODUCT CODE
■ Black	D05.110
■ Pink	D05.110P
■ Purple	D05.110PL



## TOPAZ ANEROID SPHYGMOMANOMETER

- 1 tube system
- Combined bellows & gauge tube
- Easy one-handed use
- Trigger button or screw valve pressure release
- Rotating head for left or right handed users
- 54mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult black nylon calibrated cuff
- Zip-up pouch for convenient storage
- 1 year guarantee

DESCRIPTION	PRODUCT CODE
Topaz Sphygmomanometer	D05.115



## TOPAZ DELUXE ANEROID SPHYGMOMANOMETER

- 1 tube system
- Combined bellows & gauge tube
- Easy one-handed use
- Trigger button or screw valve pressure release
- Robust metal gauge casing
- Rotating head for left or right handed users
- 58mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult black nylon calibrated cuff
- Zip-up pouch for convenient storage
- 2 year guarantee

DESCRIPTION	PRODUCT CODE
Topaz Deluxe Sphygmomanometer	D05.120



## TOPAZ MULTI-CUFF ANEROID SPHYGMOMANOMETER

- 1 tube system
- Combined bellows & gauge tube
- Easy one-handed use
- Trigger button or screw valve pressure release
- Rotating head for left or right handed users
- 54mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with large adult, adult & child cuffs
- Zip-up pouch for convenient storage
- Supplied in a presentation box
- 1 year guarantee

DESCRIPTION	PRODUCT CODE
Topaz Multi-cuff Sphygmomanometer	D05.130



## CORAL SHOCK PROOF ANEROID SPHYGMOMANOMETER

- 1 tube system
- Combined bellows & gauge tube
- Easy one-handed use
- Trigger button or screw valve pressure release
- Rotating head for left or right handed users
- Cased within a rubberised shock-proof ring
- Strong nickel-plated metal case
- 59.5mm large dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult black nylon calibrated cuff
- Zip-up pouch for convenient storage
- 2 year guarantee

DESCRIPTION	PRODUCT CODE
Coral Sphygmomanometer	D05.160





### JADE DESK/WALL ANEROID SPHYGMOMANOMETER

- 2 tube system
- Separate bellows & gauge tubes
- 152mm large dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Desk or wall mounted (supplied with wall bracket)
- Storage for cuff & tubing behind the dial
- 1 year guarantee

DESCRIPTION	PRODUCT CODE
Jade Desk/Wall Sphygmomanometer	D05.150



### OPAL ANEROID SPHYGMOMANOMETER (FLOOR MODEL)

- 2 tube system
- Separate bellows & gauge tubes
- Stand model
- 147mm large dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Heavy duty plastic feet
- Castors with brakes
- 5 wheels with 360° swivel for ease
- Adjustable height from 95 to 135cm
- Supplied with an adult black nylon calibrated cuff
- Storage for cuffs, tubes & in 1 year guarantee


DESCRIPTION	PRODUCT CODE
Opal Sphygmomanometer	D05.170



### SPHYGMOMANOMETER CUFFS

- Nylon cuffs
- Latex-free
- Push on tube connection (no adaptor required)
- Various sizes from infant to bariatric patients

DESCRIPTION	PRODUCT CODE
Obese Cuff (40.6 - 66cm), 1 Tube	D06.231
Obese Cuff (40.6 - 66cm), 2 Tubes	D06.232
Large Adult Cuff (33.3 - 51cm), 1 Tube	D06.201
Large Adult Cuff (33.3 - 51cm), 2 Tubes	D06.202
Adult Cuff (25.4 - 40.6cm), 1 Tube	D06.221
Adult Cuff (25.4 - 40.6cm), 2 Tubes	D06.222
Child Cuff (18.4 - 26.7cm), 1 Tube	D06.251
Child Cuff (18.4 - 26.7cm), 2 Tubes	D06.252
Infant Cuff (13 - 19cm), 1 Tube	D06.261
Infant Cuff (13 - 19cm), 2 Tubes	D06.262



## RUBY SINGLE HEAD STETHOSCOPE

- Quality
- Affordable
- Single head (47mm diameter)
- Ideal for use with a sphygmomanometer cuff
- Aluminium alloy chestpiece
- PVC tube
- 1 year guarantee

COLOUR	PRODUCT CODE
Blue	TD1-13-101
Green	TD1-13-104
Red	TD1-13-105
Black	TD1-13-107
Pink	TD1-13-109



## RUBY DUAL HEAD STETHOSCOPE

- Quality
- Affordable
- Dual head chestpiece (47mm diameter)
- Excellent sound transmission
- Lightweight construction
- Aluminium alloy chestpiece
- PVC tube
- 1 year guarantee

COLOUR	PRODUCT CODE
Blue	TD1-13-201
Green	TD1-13-204
Red	TD1-13-205
Black	TD1-13-207
Purple	TD1-13-208
Pink	TD1-13-209



## DUAL HEAD TEACHING STETHOSCOPE

- Dual head chestpiece (47mm diameter)
- Designed specifically for training purposes
- Two binaural headsets coupled with quality tubing
- Aluminium alloy chestpiece
- Black PVC tube
- 1 year guarantee
- Allows the teacher & student to listen to heart, chest & lung sounds simultaneously.

DESCRIPTION	PRODUCT CODE
Dual Head Teaching Stethoscope	TDI-13-217



### TWIN TUBE STETHOSCOPE (SPRAGUE RAPPAPORT)

- Dual head chestpiece
- Two tubes connecting the chestpiece to the binaural (45.5mm diameter)
- Provides excellent acoustics in each ear
- Interchangeable chestpieces for adult & paediatrics included
- Zinc alloy chestpiece
- PVC tube
- 1 year guarantee

COLOUR	PRODUCT CODE
Blue	TD1-13-302
Green	TD1-13-304
Red	TD1-13-305
Orange	TD1-13-306
Black	TD1-13-307
Purple	TD1-13-308
Pink	TD1-13-309

### DIAMOND ADULT STETHOSCOPE

- Dual head chestpiece
- A high quality 48mm diameter chestpiece with a floating diaphragm for excellent acoustics
- Contains internal binaural springs
- Plastic non-chill sleeves for patient comfort
- Supplied with a spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Stainless steel chestpiece
- PVC tube
- 3 year guarantee

COLOUR	PRODUCT CODE
Navy Blue	TD1-13-702
Burgundy	TD1-13-703
Hunter Green	TD1-13-704
Black	TD1-13-707





## DIAMOND PAEDIATRIC STETHOSCOPE

- Black
- Dual head chestpiece
- Paediatric chestpiece (40mm diameter)
- Floating diaphragm for excellent acoustics
- Contains internal binaural springs
- Plastic non-chill sleeves for patient comfort
- Supplied with spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Stainless steel chestpiece
- PVC tube
- 3 year guarantee

DESCRIPTION	PRODUCT CODE
Diamond Paediatric Stethoscope	TD1-13-717



## DIAMOND INFANT STETHOSCOPE

- Black
- Dual head chestpiece
- Infant chestpiece (30mm diameter)
- Floating diaphragm for excellent acoustics
- Contains internal binaural springs
- Plastic non-chill sleeves for patient comfort
- Supplied with spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Stainless steel chestpiece
- PVC tube
- 3 year guarantee

DESCRIPTION	PRODUCT CODE
Diamond Infant Stethoscope	TD1-13-727



## DIAMOND NEONATAL STETHOSCOPE

- Black
- Dual head chestpiece
- Neonatal chestpiece (24mm diameter)
- Floating diaphragm for excellent acoustics
- Contains internal binaural springs
- Plastic non-chill sleeves for patient comfort
- Supplied with a spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Stainless steel chestpiece
- 3 year guarantee

DESCRIPTION	PRODUCT CODE
Diamond Neonatal Stethoscope	TD1-13-737

### CARDIOLOGY ADULT STETHOSCOPE



- Dual head chestpiece
- 48mm diameter chestpiece
- Deep cone-shaped bell for excellent acoustic sensitivity
- High quality tube for resistance to ambient interference
- A-shaped tube design features a binaural with reinforced stainless steel inner-spring for durability
- Features an anti-chill rim for patient comfort
- Supplied with a spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Stainless steel chestpiece
- PVC tube
- 3 year guarantee

COLOUR	PRODUCT CODE
■ Navy Blue	TD1-13-902
■ Black	TD1-13-907

### DELUXE CARDIOLOGY ADULT STETHOSCOPE



- Ergonomic & unique design
- 50mm diameter chestpiece
- Quality tube for excellent acoustic performance
- Comfortable & flexible binaural
- Durable, soft silicone ear tips
- Supplied with a spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Zinc alloy chestpiece
- PVC tube
- 3 year guarantee

COLOUR	PRODUCT CODE
■ Black	TD1-13-807

### PINARD STETHOSCOPE



- For listening to fetal heart sounds
- Available in plastic, aluminium or wood
- Approximately 14.5cm

COLOUR & SIZE	PRODUCT CODE
■ Plastic Black	03.350.10
■ Plastic Red	03.350.15
■ Plastic Pink	03.350.20
■ Aluminium	03.360.20
■ Wood	03.370.20





## EMERALD SPHYGMOMANOMETER & SPRAGGUE STETHOSCOPE SET

### Sphygmomanometer

- Separate bellows & gauge tubes
- Spring loaded clip-on rear of dial for one-handed use
- 48mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult nylon calibrated cuff

### Stethoscope

- Dual head chestpiece
- 2 tubes connecting the chestpiece to the binaural
- Interchangeable chestpieces for adult & paediatrics incl.

DESCRIPTION	PRODUCT CODE
Black Emerald Aneroid	
■ Sphygmomanometer & Black Twin Tube Sprague Rappaport Stethoscope Set	D05.305
Pink Emerald Aneroid	
■ Sphygmomanometer & Pink Twin Tube Sprague Rappaport Stethoscope Set	D05.305P
Purple Emerald Aneroid	
■ Sphygmomanometer & Purple Twin Tube Sprague Rappaport Stethoscope Set	D05.305PL



## MIDWIFERY SET PINK (COMPLETE)

### Includes:

- Stethoscope
- Sphygmomanometer
- Pinard
- Tourniquet
- Fob watch
- Pen torch
- Face shield
- Dressing scissors
- Artery forceps
- Hand sanitiser
- Utility scissors
- Obstetric wheel

DESCRIPTION	PRODUCT CODE
Midwifery Set Pink (Complete)	D05.905
Pink Plastic Case	D05.950



## NURSES SET PINK (COMPLETE)

### Includes:

- Stethoscope
- Sphygmomanometer
- Tourniquet
- Fob watch
- Pen torch
- Face shield
- Dressing scissors
- Artery forceps
- Hand sanitiser
- Utility scissors

DESCRIPTION	PRODUCT CODE
Nurses Set Pink (Complete)	D05.910
Pink Plastic Case	D05.950



### DISPOSABLE PEN TORCH



- Single-use
- Pocket clip
- Plastic sealed casing
- Non-replaceable battery
- Features a pupil gauge
- Integrated pocket clip, when pressed powers the torch
- Supplied in packs of six

DESCRIPTION	PRODUCT CODE
Disposable pen torch	D85.100



### PLASTIC REUSABLE PEN TORCH

- Reusable (batteries can be replaced)
- Plastic casing
- Integrated pocket clip, when pressed powers the torch
- Supplied with 2\*AAA batteries
- Provided in a blister presentation pack

COLOUR	PRODUCT CODE
■ Blue	D85.005
■ Green	D85.010
■ Red	D85.015
□ White	D85.020
■ Pink	D85.025
■ Purple	D85.030



### DELUXE STAINLESS STEEL REUSABLE PEN TORCH

- Reusable (batteries can be replaced)
- Stainless steel casing
- Features a pupil gauge
- Pocket clip
- Supplied with 2\*AAA batteries
- Provided in a blister presentation pack

DESCRIPTION	PRODUCT CODE
Deluxe pen torch	D85.300



## SILICONE FOB WATCH



Pin clasp



Battery save clip

- New improved design
- Japanese Quartz movement
- Pin clasp closure
- Made from silicone based material
- Ideal for infection control
- Contains 1 x watch battery (SR626SW)
- 2 year battery guarantee

COLOUR	PRODUCT CODE
<span style="color: pink;">■</span> Pink	TDM-FW-500
<span style="color: blue;">■</span> Blue	TDM-FW-505
<span style="color: green;">■</span> Green	TDM-FW-510
<span style="color: red;">■</span> Red	TDM-FW-515
<span style="color: white;">□</span> White	TDM-FW-520
<span style="color: black;">■</span> Black	TDM-FW-525
<span style="color: purple;">■</span> Purple	TDM-FW-530
<span style="color: black;">■</span> Pack of 14 (2 of each colour)	TDM-FW-545



## CARABINER WATCH

- Carabiner clip-on style
- Robust metal casing
- Clip to prevent battery engagement prior to use
- Japanese Quartz movement
- Contains 1 x watch battery (SR626SW)
- 2 year battery guarantee

DESCRIPTION	PRODUCT CODE
Carabiner fob watch	TDM-FW-550



## UTILITY SCISSORS

- Tuff cuts
- Stainless steel blades
- Tough/durable plastic handle
- Available in a variety of colours
- Lower blade has a serrated edge & is hooked & rounded at the lip for ease of use
- Available in two sizes, 19cm/7.5" & 15cm/6"
- Pack of 10

COLOUR & DESCRIPTION	PRODUCT CODE
Utility Scissors 7.5" Green	12.711.06
Utility Scissors 7.5" Yellow	12.711.11
Utility Scissors 7.5" Orange	12.711.16
Utility Scissors 7.5" Blue	12.711.21
Utility Scissors 7.5" Black	12.711.31
Utility Scissors 7.5" Red	12.711.36
Utility Scissors 7.5" Pink	12.711.41
Utility Scissors 7.5" Purple	12.711.46

COLOUR & DESCRIPTION	PRODUCT CODE
Utility Scissors 6" Green	12.751.06
Utility Scissors 6" Blue	12.751.21
Utility Scissors 6" Black	12.751.31
Utility Scissors 6" Red	12.751.36
Utility Scissors 6" Pink	12.751.41
Utility Scissors 6" Purple	12.751.46



## TOURNIBUTTON TOURNIQUET



- Single-use to reduce the risk of cross infection
- Supplied in a convenient dispensing box
- Easy to reposition & tighten with a button/stud fastener
- Adult & Paediatric use
- 2 bands can be joined for bariatric use
- Material: TPE band
- Latex-free
- Band size: 38.5cm x 2.3cm
- Box of 100

DESCRIPTION	PRODUCT CODE
TourniButton Single-Use Tourniquets (Button Type)	D80.200



## TOURNIBAND TOURNIQUET



- Single-use to reduce the risk of cross infection
- Supplied on a perforated roll
- Stretch band
- Supplied in dispenser pack
- Latex-free
- Material: TPE
- Band size: 45cm x 2.5cm (approx.)
- Available in boxes of 25 or 100

DESCRIPTION	PRODUCT CODE
TourniBand Single Use Tourniquets Box 25	D80.100
TourniBand Single Use Tourniquets Box 100	D80.101



## TOURNICLIP QUICK RELEASE TOURNIQUET

- Reusable
- Non pinch for patient comfort
- Quick clip application & button release for ease of use
- Available in a variety of colours
- Latex-free

COLOUR	PRODUCT CODE
Blue	D80.005
Green	D80.010
Red	D80.015
Pink	D80.020
Purple	D80.025
Yellow	D80.030



## EVAC-S SINGLE-USE SUCTION DEVICE



- Manual handheld suction device
- Completely single-use
- Easy to assemble
- Product is ready to use after connecting the catheter
- After use the device can be placed into the polybag, resealed & the product label can be peeled to reveal a biohazard warning label
- Vacuum (max): 500mmHg
- Peak free airflow: >20ml
  - Disposable container volume: 300ml
  - External diameter of cap (suction port): 17mm
  - External diameter of catheters: 9mm Adult, 10fg Paediatric
- Operating environmental temperature: -20°C to +50°C
- Storage environmental temperature: -40°C to +60°C
- Net weight: 270g



The box contents are:

- Suction pump handle
- 300ml canister
- Adult catheter
- Paediatric catheter
- Self-seal disposal polybag
- Bio-hazard label ready to be displayed for disposal
- Multilingual product instructions
- Multilingual product labeling
- Supplied in boxes of 5

DESCRIPTION	PRODUCT CODE
EVAC-S suction device	TDM-ESD-100

## CPR RESCU-MASK



- Single-use
- Replacement filters available for training purposes

Supplied with:

- One-way bacterial & viral filter
- Built-in O2 port
- Elasticated head strap
- Hard plastic carry case

DESCRIPTION	PRODUCT CODE
CPR Rescu-Mask	TDM-PM-250
One-way valve / filter (x5)	TDM-PM-110



## FACE SHIELD IN KEYRING POUCH



- Single-use
- Available in various coloured keyring pouches
- Convenient & effective protection from cross-infection during mouth-to-mouth resuscitation
- One-way valve filter

COLOUR	PRODUCT CODE
Blue	TDM-PM-051
Black	TDM-PM-052
Purple	TDM-PM-053
Red	TDM-PM-054
Yellow	TDM-PM-056
Orange	TDM-PM-057
Green	TDM-PM-058
Pink	TDM-PM-059
Pack of 40 (5 of each colour)	TDM-PM-045



## TIMESCO DISPOSABLE MANIKIN FACE SHIELDS



- Disposable manikin face shields for use in mouth-to-mouth resuscitation/CPR training situations
- Helps ensure student safety
- Supplied in rolls of 50

DESCRIPTION	PRODUCT CODE
Timesco Disposable Manikin Face Shields Roll of 50	TDM-PM-175



## FACE SHIELD IN SACHET



- Face shield for personal protection during mouth-to-mouth resuscitation
- Supplied in easy tear sachet
- Paper filter

DESCRIPTION	PRODUCT CODE
Face shield in sachet	TDM-PM-165

## GO GREEN.

We accept that our daily activities may have both a positive & negative impact on the environment. We are committed to managing our operations, products & services in such a way to preserve, protect & enhance the environment. Our goal is to minimise this effect & to reduce our carbon footprint wherever possible.

We aim to achieve this by:

- Ensuring that every aspect of our activities is conducted in accordance with environmental policies.
- Minimising the consumption of natural resources & energy, whilst consuming material goods in moderation.
- Reducing the creation of waste by improving practices & recycling materials whenever practical.
- Continually seeking to improve environmental practices & performance & develop more environmentally friendly products.
- Informing & educating all Timesco staff about environmental issues in the context of the business through appropriate training & information sharing.

Timesco Healthcare Ltd., 3 Carnival Park, Carnival Close,  
Basildon, Essex, SS14 3WN, United Kingdom  
T: +44 (0)1268 297 700 E: [info@timesco.com](mailto:info@timesco.com) W: [timesco.com](https://www.timesco.com)



**TIMESCO**