

callisto
TIMESCO™

The Timesco Callisto brand is trusted and respected due to its high level of clinical performance.

An extensive range of single use blades including specialist blades, all with exceptional specifications including; a rigid metal spatula, secure fitting and no touch hook to reduce the risk of cross infection.



CE marked

**ISO
7376**

ISO 7376
compatible



Single Use



Supplied with
multilingual
instructions



Latex Free



Compatible
handles are
colour coded

SPECIFICATION



FEATURES



Size 3

MACINTOSH

Product Code	Description	Measurements*
DS.3940.150.05	Callisto Macintosh Size 0	80mm x 13mm
DS.3940.150.10	Callisto Macintosh Size 1	90mm x 18mm
DS.3940.150.15	Callisto Macintosh Size 2	100mm x 18mm
DS.3940.150.20	Callisto Macintosh Size 3	130mm x 18mm
DS.3940.150.25	Callisto Macintosh Size 4	153mm x 18mm



Size 1

MILLER

Product Code	Description	Measurements*
DS.3940.185.03	Callisto Miller Size 00	65mm x 12mm
DS.3940.185.05	Callisto Miller Size 0	79mm x 12mm
DS.3940.185.10	Callisto Miller Size 1	105mm x 13mm
DS.3940.185.15	Callisto Miller Size 2	158mm x 15mm
DS.3940.185.20	Callisto Miller Size 3	197mm x 16mm
DS.3940.185.25	Callisto Miller Size 4	208mm x 18mm



Size 1

ROBERTSHAW

Product Code	Description	Measurements*
DS.3940.235.05	Callisto Robertshaw Size 0	94mm x 10mm
DS.3940.235.10	Callisto Robertshaw Size 1	112mm x 12mm



Size 1

SEWARD

Product Code	Description	Measurements*
DS.3940.240.10	Callisto Seward Size 1	107mm x 13mm

Timesco Healthcare Ltd Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800
E info@timesco.com | timesco.com

L.L.054/1-W E & O.E Timesco Dec 2017.

*Length x Height. Height is measured from the top of the block to the top of the blade. Measurements rounded to the nearest millimeter.