

optima^{CLX}

TIMESCO™

Optima CLX are reusable blades with a removable glass fibre light guide.

The light guide easily clicks in and out without the use of tools to allow cleaning and replacement of the core, thus reducing the cost of ownership and maintenance.



CE marked



Fully autoclaveable



Latex Free



Extensive guarantees

**ISO
7376**

ISO 7376 compatible



Supplied with multilingual instructions



Compatible handles are colour coded

SPECIFICATION

FEATURES

The fibre core is positioned to offer the least obstructive view.

Highly polished 4mm replaceable, glass fibre core for focused illumination. It has a one piece snap in/out fibre core for easy cleaning and replacement. No tools required.



Clear product information etched on the blade.

3 metal ball bearings for secure fitting.



A green, fibre optic ISO 7376 fitting metal block.

Atraumatic distal end.

Fully autoclaveable metal blade with an anti-glare satin finish. Covered by a 5 year guarantee or 4,000 autoclaves, whichever is sooner. See website for full terms and conditions.



Size 3

MACINTOSH

Product Code	Description	Measurements*
2970.150.05	Optima CLX Macintosh Size 0	79mm x 15.5mm
2970.150.10	Optima CLX Macintosh Size 1	92mm x 22mm
2970.150.15	Optima CLX Macintosh Size 2	100mm x 21mm
2970.150.20	Optima CLX Macintosh Size 3	131mm x 23.5mm
2970.150.25	Optima CLX Macintosh Size 4	155mm x 25mm



Size 1

MILLER

Product Code	Description	Measurements*
2970.185.03	Optima CLX Miller Size 00	66.5mm x 11.5mm
2970.185.05	Optima CLX Miller Size 0	77mm x 11.5mm
2970.185.10	Optima CLX Miller Size 1	101.5mm x 12mm
2970.185.15	Optima CLX Miller Size 2	153mm x 15mm
2970.185.20	Optima CLX Miller Size 3	197mm x 15mm
2970.185.25	Optima CLX Miller Size 4	205mm x 17mm

In addition to the Optima CLX blades above, we also supply speciality blades, including Robertshaw, Seward and Wisconsin patterns.

Spare Optima CLX cores are also available.

Please call 01268 297 710 or visit www.timesco.com for further details.

Timesco Healthcare Ltd Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800
E info@timesco.com | timesco.com

L.L.059/1-W E & O.E Timesco Feb 2018.

*Length x Height. Height is measured from the top of the block to the top of the blade. Measurements rounded to the nearest 0.5 millimeter.