

siriusXL
TIMESCO™

Sirius XL are reusable blades with an enclosed, non-replaceable fixed glass fibre light guide. The system requires no disassembly for cleaning.



CE marked



Fully autoclaveable



Latex Free



Extensive guarantees

**ISO
7376**

ISO 7376 compatible



Supplied with multilingual instructions



Compatible handles are colour coded

SPECIFICATION

FEATURES

The non-replaceable glass fibre bundle is fully enclosed in the blade which allows easy cleaning.

The fibre core is positioned to offer the least obstructed view.

Atraumatic distal end.

Fully autoclaveable metal blade, covered by a 5 year guarantee or 4,000 autoclaves, whichever is sooner. See website for full terms and conditions.

Green ISO 7376 fitting metal block compatible with fibre optic handle.

Clear product information etched on the blade.

3 metal ball bearings for secure fitting.





Size 3

MACINTOSH

Product Code	Description	Measurements*
2955.150.05	Sirius XL Macintosh Size 0	83mm x 20mm
2955.150.10	Sirius XL Macintosh Size 1	95mm x 20mm
2955.150.15	Sirius XL Macintosh Size 2	115mm x 24mm
2955.150.20	Sirius XL Macintosh Size 3	135mm x 24mm
2955.150.25	Sirius XL Macintosh Size 4	160mm x 24mm



Size 1

MILLER

Product Code	Description	Measurements*
2955.185.03	Sirius XL Miller Size 00	65mm x 18mm
2955.185.05	Sirius XL Miller Size 0	77mm x 18mm
2955.185.10	Sirius XL Miller Size 1	102mm x 18mm
2955.185.15	Sirius XL Miller Size 2	150mm x 19mm
2955.185.20	Sirius XL Miller Size 3	195mm x 19mm
2955.185.25	Sirius XL Miller Size 4	207mm x 22mm

Timesco Healthcare Ltd Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800
E info@timesco.com | timesco.com

L.L.057/1 E & O.E Timesco Feb 2018.

*Length x Height. Height is measured from the top of the block to the top of the blade. Measurements rounded to the nearest millimeter.