

The compact, portable and water-resistant Pentax Airway Scope video laryngoscope provides safe, fast and effective tracheal intubation. High resolution colour images and a unique targeting system facilitate precise intubation from virtually any position to minimise patient trauma and ensure accurate placement of endotracheal tubes using compatible PBlades.



CE Marked



Supplied with
multilingual
instructions



LCD colour screen



LED light

SPECIFICATION

Up to a 30 mega-pixel high definition clear colour LCD screen provides vivid images of the larynx and glottis.

A cross-hair display on the screen supports accurate and reliable intubation.

This waterproof model can be used outdoors even in wet conditions.



The Pentax Airway Scope (AWS-S200) weighs just 235g making it the lightest model yet.

Benefiting from low power consumption and designed for up to 1 hour continuous use on AA batteries.

FEATURES



The target symbol on the monitor facilitates safer intubation, even for less experienced users.

Connect to medical monitors or video recording devices to simultaneously show footage from the scope on a larger monitor or to record.

With the variety of PBLADE family, securing the airway through tracheal intubation and suction of the oral cavity can be performed smoothly without tracheal damage.



PENTAX AIRWAY SCOPE

Product Code	Description
6000.100.00	Pentax Airway Scope

Equipped with a tube channel so no stylett (guide wire) required.

Sterile and single use to reduce the risk of cross infection.



The tube channel allows easy intubation even after insertion of a tube into the oral cavity.

Compatible with tracheal tubes of an outer diameter of 8.5mm to 11.

Hole to accommodate suction catheters 4.0mm (12F) or less for suction of the oral cavity.



Neonate PBLADE



Paediatric PBLADE



Thin PBLADE



Standard PBLADE

The unique design of the rigid, sterile, single use PBlade used with the Pentax Airway Scope provides an unparalleled clear view for safe, smooth & accurate intubation whilst eliminating the risk of cross infection.

CE CE marked



Supplied with
multilingual
instructions



Single Use

STERILE EO

Sterilised

PENTAX PBLADES

Product Code	Description	Measurements
6000.110.20	PBLADE (ITL-NL NEONATE)	H148 x W33 x D53mm
6000.110.15	PBLADE (M-ITL-PL PAEDIATRIC)	H140 x W33 x D65mm
6000.110.10	PBLADE (M-ITL-TL THINNER)	H130 x W44 x D90mm
6000.110.05	PBLADE (MITL-SL STANDARD)	H131 x W49 x D96mm

Coming soon

Double Lumen PBlade. Contact
sales@timesco.com for further
information.



Neonate PBLADE



Paediatric PBLADE



Thin PBLADE



Standard PBLADE

Timesco Healthcare Ltd Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | **F** +44 (0) 1268 297 800
E info@timesco.com | timesco.com

L.L.063/1 E & O.E Timesco April 2018