

REUSABLE FACE MASKS

TIMESCO™

The Timesco reusable silicone face masks are transparent for optimum patient visibility and cushioned providing exceptional patient comfort. They are fully autoclavable and available in a full range of sizes from premature to large adult.



CE Marked

ISO

ISO 10651
ISO 5356-1
ISO 10993



Fully
autoclavable



Supplied with
multilingual
instructions



Latex Free

SPECIFICATION

Transparent face piece allows visual patient monitoring at all times and visual checks for bleeding and vomiting.

Connectors are compatible with all standard catheter mounts and breathing systems.



The moulded texture of the rim ensures maximum patient comfort and a tight seal.

Available in premature, neonate, infant, paediatric, medium, standard and large adult sizes.

FEATURES



TDM-MP-1500



TDM-MP-1513



TDM-MP-1515

Product Code	Description	Size	NPC Code
TDM-MP-1506	Premature	Size 00	FDD5161
TDM-MP-1500	Neonate	Size 0	FDD668
TDM-MP-1501	Infant	Size 1	FDD669
TDM-MP-1502	Paediatric	Size 2	FDD670
TDM-MP-1513	Small Adult	Size 3	FDD671
TDM-MP-1514	Medium Adult	Size 4	FDD1055
TDM-MP-1515	Large Adult	Size 5	FDD1056

Timesco Healthcare Ltd Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800
E info@timesco.com | timesco.com

L.L.066/1 E & O.E Timesco January 2019