



The Timesco new & improved XLED handles provide exceptionally bright white illumination from a replaceable LED bulb supplying consistent, flicker free light & maximising battery life due to low power consumption. They are available in medium, small & stubby size and are compatible with reusable & single use blades.



CE marked



Fully autoclavable (excluding bulb)



Supplied with multilingual instructions



Extensive guarantee



Latex-free



ISO 7376 compatible



Compatible blades are colour coded



Replaceable, long life LED bulb

SPECIFICATION

A polished metal collar reflects and directs light up towards the fibre core of the blade for maximum illumination.

Colour coded band to clearly identify it is ISO 7376 fitting compatible with green fitting blades and clearly marked with the XLED logo, product code, CE mark and LOT number.

The metal handle is fully autoclavable (excluding bulb) and covered by a 5 year guarantee or 4,000 autoclaves, whichever is sooner. Refer to our website www.timesco.com for full terms and conditions.

Available in medium, small & stubby size and are compatible with green fitting reusable & single use blades.



High performance & reliability from a super-bright replaceable 3v LED bulb providing excellent, flicker free, white illumination.

The bulb provides 3 x more lux output than other leading brands with an exceptional lux output of 40,000 (at bulb) and a consistent colour temperature of 5,000(k).

The XLED handle works with dry cell or rechargeable batteries without requiring a voltage of 3.5v. Intelligent electronics control the light intensity to provide consistent, reliable & flicker free illumination, even at low battery voltages.

Knurled handle design provides excellent grip. Secure screw cap for easy access to batteries.

FEATURES



Small



Medium



Stubby

Product Code	Description	Measurements*	NPC Code
3000.500.05	XLED Handle, Small	146mm x 19mm	FSM718
3000.500.10	XLED Handle, Medium	146mm x 29mm	FSM717
3000.508.10	XLED Handle, Stubby	111mm x 32mm	FSM719
3000.950.20	Spare Bulbs (Pack of 6)	-----	-----

Timesco Healthcare Ltd Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800

E info@timesco.com | timesco.com

L.L.050/4 E & O.E Timesco August 2019.

*Height x Diameter. Measurements rounded to the nearest millimeter.