

callisto
TIMESCO™

Callisto is the U.K.'s No.1 single use laryngoscope system from Timesco. Our Callisto sets benefit from exceptional specifications in a convenient ready to use set and are available in an extensive range.

They offer a true low cost alternative to reusable systems, whilst reducing the risk of cross infection and eliminating maintenance costs.



CE marked



Single-use



Latex-free

**ISO
7376**

ISO 7376
compatible



Supplied with
multilingual
instructions



Fibre-optic,
green fitting

SPECIFICATION

Rigid metal spatula looks and feels like a reusable fibre optic blade.

Highly polished 4mm acrylic core for maximum illumination.

Available in a comprehensive range of convenient ready to use sets.

Preloaded with alkaline button cell batteries, which are quick and easy to remove for safe battery disposal.



Metal hook bar offering exceptional strength for a single use handle and contains 2 metal ball bearings for secure fitting.

Powerful bright LED bulb for excellent illumination with an average minimum output of 11,000 lux when measured 2cm from the bulb.

Features a low volume plastic construction to reduce plastic waste and cost of disposal.

FEATURES



CALLISTO MACINTOSH BLADE & LED PRELOADED COMPACT HANDLE SET

Product Code	Description	Quantity	NPC Code
DS.3950.952.05	Callisto Blade Macintosh 0 & LED Preloaded Compact Handle	Box of 10	FSM3925
DS.3950.952.10	Callisto Blade Macintosh 1 & LED Preloaded Compact Handle	Box of 10	FSM3813
DS.3950.952.15	Callisto Blade Macintosh 2 & LED Preloaded Compact Handle	Box of 10	FSM3229
DS.3950.952.20	Callisto Blade Macintosh 3 & LED Preloaded Compact Handle	Box of 10	FSM3230
DS.3950.952.25	Callisto Blade Macintosh 4 & LED Preloaded Compact Handle	Box of 10	FSM3231



CALLISTO MILLER BLADE & LED PRELOADED COMPACT HANDLE SET

Product Code	Description	Quantity	NPC Code
DS.3950.952.28	Callisto Blade Miller 00 & LED Preloaded Compact Handle	Box of 10	FSM3814
DS.3950.952.30	Callisto Blade Miller 0 & LED Preloaded Compact Handle	Box of 10	FSM3232
DS.3950.952.35	Callisto Blade Miller 1 & LED Preloaded Compact Handle	Box of 10	FSM3233
DS.3950.952.40	Callisto Blade Miller 2 & LED Preloaded Compact Handle	Box of 10	FSM3815
DS.3950.952.45	Callisto Blade Miller 3 & LED Preloaded Compact Handle	Box of 10	FSM3816



CALLISTO ROBERTSHAW BLADE & LED PRELOADED COMPACT HANDLE SET

Product Code	Description	Quantity	NPC Code
DS.3950.952.55	Callisto Blade Robertshaw 0 & LED Preloaded Compact Handle	Box of 10	FSM3926
DS.3950.952.60	Callisto Blade Robertshaw 1 & LED Preloaded Compact Handle	Box of 10	FSM3927

Timesco Healthcare Ltd Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800
E info@timesco.com | timesco.com

L.L.073/1 E & O.E Timesco August 2019. Images shown are for illustration purposes only.