

• NEW • FETAL DOPPLER • LOW COST • HIGH PERFORMANCE •

Our new non-invasive, compact hand held fetal doppler is intended for general purpose screening of the fetal heart throughout all stages of pregnancy.

Auto power off

3.0MHz ultrasound probe

Hand held device

Portable & compact

Volume control


Non-invasive

Clear output from speaker

Detachable probe



Clear TFT screen displaying:



• AVAILABLE DIRECT FROM TIMESCO FOR ONLY £95 •

**TIMESCO**  
INSTRUMENTS FOR LIFE

# Fetal Doppler Specification Comparison Sheet



Product name & Timesco code	ChoiceM Med MD800 - TDM-MD800	Huntleigh SR2 - TDM-HT-SR2	Ultra-Tec DWR - TDM-PD1DWR
Warranty	1 Year	3 Years	3 Years
High Resolution or LCD Screen	High Resolution (TFT)	High Resolution	LCD
Coloured Screen	Yes	Yes	No
Fetal Heart Rate Display (FHR)	Yes	Yes	Yes
Fetal Heart Rate Waveform Display	Yes	Yes	No
Pulse Indicator	Yes	Yes	Yes
Probe	High Sensitivity, 3.0MHz	High Sensitivity, 2.0MHz	High Sensitivity, 3.0MHz
Measuring Range	50bpm - 210bpm	50bpm - 210bpm	50bpm - 210bpm
Accuracy	+/- 2bpm	+/- 3bpm	+/- 1bpm
Resolution Ratio	1bpm	1bpm	1bpm
Dimensions (L x W x H)	130mm x 90mm x 38mm	140mm x 75mm x 30mm	150mm x 75mm x 31mm
Weight (without batteries)	190g	310g	295g
Waterproof	No	Yes (Probe Only)	Yes (Probe Only)
Auto Power Off	Yes (1 minute)	Yes (1 minute)	Yes (5 minutes)
Built In Speaker	Yes (Adjustable)	Yes	Yes (Adjustable)
Data Storage	No	Yes	No
Software Compatible	No	Yes	No
Safety Standards	IEC60601-1:2006 & CE	IEC60601-1:2006 & CE	IEC60601-1:2006 & CE
Supplied With Carry Case	No	Yes	Yes
Supplied With Gel	No	Yes	Yes
Supplied With Batteries	Yes (2 x AA Batteries)	Yes (2 x AA Batteries)	Yes (1 x 9V Battery)
Low Battery Indicator	Yes	Yes	Yes
Battery Life	8 Hours Continuous Use	8 Hours Continuous Use	6 Hours Continuous Use
Price	£95.00 (ex. VAT)	£447.00 (ex. VAT)	£323.00 (ex. VAT)