


TIMESCO'S RANGE OF THERMOMETERS



AVAILABLE VIA NHS SUPPLY CHAIN



Available via NHS Supply Chain Framework Agreement



1 FWH225

- Measurement time 60 seconds
- 1 year warranty
- Rigid tip
- Auto shut off
- Celsius & Fahrenheit readings
- Use with [FWH226](#) Probe Covers



2 FWH224 / FWH123

- Measurement time 10 seconds
- 1 year warranty
- Rigid tip (224) / Flexible tip (123)
- Auto shut off
- Celsius & Fahrenheit readings
- Use with [FWH226](#) Probe Covers



3 FWH120

- Measurement time 1 second
- 1 year warranty
- Auto shut off
- 2 x AAA batteries
- Celsius readings
- Use with [FWH121](#) Probe Covers



4 FWH122

- 1 year warranty
- Non contact
- Non invasive
- Infection control friendly
- No spare ear probes required
- No ongoing costs

	NPC Code	Product Code	Description	Price (inc VAT)
1	FWH225	D60.215	Eco Digital Thermometer Rigid	£1.41
2	FWH224	D60.220	Rapid Digital Thermometer Rigid Tip (Pack 10)	£17.14
2	FWH123	D60.225	Rapid Digital Thermometer Flexible Tip (Pack 10)	£20.57
3	FWH120	D60.200	Digital Ear Thermometer	£17.43
4	FWH122	D60.300	Thermofinder Non-Contact Infra-Red Thermometer	£44.15
	FWH226	D60.230	Digital Thermometer Covers (sheaths) (Box 100)	£1.60
	FWH121	D60.210/10	Ear Thermometer Probe Covers (Pack 200)	£1.41

Available via NHS Supply Chain: www.supplychain.nhs.uk or contact Timesco on +44 (0) 1268 297 710

Timesco Healthcare Ltd 3 Carnival Park Basildon Essex SS14 3WN England

T +44 (0) 1268 297 710 | F +44 (0) 1268 297 801 | E sales@timesco.com | timesco.com

L.GM.009/2 E & O.E Timesco June 2018.

TIMESCO
INSTRUMENTS FOR LIFE