

# Timesco Product Guarantee

## Reusable Laryngoscopes and Surgical Instruments

The respective guarantee will be the number of years or the number of autoclave cycles, whichever is sooner and is from the date of purchase. Autoclave cycles are approved between 134°C and 137°C for 3 minutes holding time.

If any of Timesco's reusable products crack, bend, deform, corrode or fail whilst being used for the intended purpose within the guarantee period, we will repair it or replace it free of charge.

Range	Guarantee	Number of autoclave cycles
<b>Reusable Laryngoscope Blades and Handles:</b>		
Sirius, Optima CLX, Eclipse and Orion laryngoscope blades <i>(excluding bulbs)</i>	5 years	4,000*
XLED, Xenon and Orion laryngoscope handles <i>(excluding bulbs)</i>	5 years	4,000*
Orion bulbs	6 months	100
<b>Reusable Surgical Instruments:</b>		
Reusable surgical instruments (excluding the sections of the products listed below)	5 years	4,000
• The cutting, shearing, grinding and rasping edges (excluding Tungsten Carbide)	1 year	800
• The cutting, shearing, grinding and rasping edges (Tungsten Carbide)	2 years	1,600
Reusable Electrosurgery cables	1 year	100
Reusable Bipolar and Monopolar forceps	1 year	40

\*All blade and handle combinations will exceed the minimum requirement of 500 LUX after 4,000 autoclave cycles in accordance with the ISO 7376 standard.

For more information and the terms and conditions please visit [www.timesco.com](http://www.timesco.com) or call 01268 297 700.

Timesco Healthcare Ltd., 3 Carnival Park, Basildon, Essex, SS14 3WN England.  
 T: +44 (0)1268 297710 | F: +44 (0)1268 297800  
 E: [sales@timesco.com](mailto:sales@timesco.com) | [www.timesco.com](http://www.timesco.com)

