

RESUSCITATOR (BVM) OVERVIEW

Timesco's Resuscitators (BVMs) are single patient use and supplied with oxygen tubing, oxygen reservoir bag and suitably sized face masks.

The resuscitators (BVMs) are supplied concertinaed for space saving and pre-assembled for immediate use with glued components (except the interchangeable mask).

They are available in a range of sizes:
Adult - 1000ml or 1500ml | Child - 500ml | Infant - 280ml

Latex free and do not contain DEHP.

Product Code	Description	NHS Supply Chain Code
TDM-RES-0280	Infant 280ml BVM - Size 0 Mask	FDE1627
TDM-RES-0550	Child 550ml BVM - Size 1 & 2 Masks	FDE1626
TDM-RES-1000	Adult 1000ml BVM - Size 5 Mask	FDE1628
TDM-RES-1500	Adult 1500ml BVM - Size 5 Mask	FDE1629
TDM-RES-1500W	Adult 1500ml BVM with Pop-off Valve - Size 5 Mask	FDE1625

Adult is available with or without 'Pop-Off' pressure relief valve (40-60cm H₂O)

Child is supplied with 35-50cm H₂O

Infant is supplied with 35-45cm H₂O



Textured grip for secure handling

Fast recoil to allow for rapid and consistent ventilation



Swivel between valve and mask which allows 360° positioning in relation to the patient

Kink resistant 2m oxygen tubing, featuring an internal star shape

Adjustable air cushion face mask for maximum patient comfort (no hook rings)

Transparent bag and mask allows complete internal visibility

QUESTIONS?

Simply get in touch via the contact details below.

T 01268 297 710 | E sales@timesco.com | Live chat: timesco.com