



# FOETAL DOPPLERS

LOW COST • WATERPROOF PROBES • GUARANTEED

★ NEW PRODUCT



★ NEW PRODUCT



**HI-BEBE (BT-200L) FOETAL DOPPLER**

- Small, portable and easy to operate
- Lightweight
- Foetal heart rate processing
- LCD display
- Built in loudspeaker
- 2 MHZ frequency
- Waterproof probe
- FHR and sound management
- Supplied with carry case and gel
- 1 year guarantee

PRODUCT CODE	DESCRIPTION
TDM-BT-200L	Hi-Bebe (BT-200L) Foetal Doppler

NPC CODE	NHS SC BAND PRICE 1	NHS SC BAND PRICE 2
FGR1023	£79.95 (1-9 units)	£71.95 (10+ units)

**FEATURES**

**WATERPROOF PROBE**

Suitable for water birth delivery. (IPX7)

**2MHZ FREQUENCY**

Supplied with a 2MHz probe.

**SUITABLE FOR ALL STAGES OF PREGNANCY**

Intended for general purpose screening of the foetal heart rate (FHR) at all stages of pregnancy.

**HI-BEBE (BT-220C) FOETAL DOPPLER**

- Small, portable and easy to operate
- Lightweight
- Foetal heart rate processing
- 2.4" LCD display with FHR trend display
- Built in loudspeaker
- 2 MHZ frequency
- Waterproof probe
- FHR and sound management
- Acoustic stimulator function
- BMI and body fat analysis
- Supplied with carry case and gel
- 1 year guarantee

PRODUCT CODE	DESCRIPTION
TDM-BT-220C	Hi-Bebe (BT-220C) Foetal Doppler

NPC CODE	NHS SC BAND PRICE 1	NHS SC BAND PRICE 2
FGR1024	£144.95 (1-9 units)	£130.50 (10+ units)

**FEATURES**

**WATERPROOF PROBE**

Suitable for water birth delivery. (IPX7)

**ACOUSTIC STIMULATOR**

Stimulate sleeping foetus to acquire more accurate FHR. High (100Hz) / Low (75Hz)

**SUITABLE FOR ALL STAGES OF PREGNANCY**

Intended for general purpose screening of the foetal heart rate (FHR) at all stages of pregnancy.

**INTERCHANGEABLE PROBE OPTION**

Supplied with a 2MHz probe.

**TREND DISPLAY**

Large colour LCD enables trend of FHR on screen.

## FOETAL DOPPLER SPECIFICATION COMPARISON SHEET



Product name / Timesco code	Hi-Bebe BT-200L / TDM-BT-200L	Hi-Bebe BT-220C / TDM-BT-220C
Guarantee	1 Year	1 Year
High resolution or LCD screen	Mono high resolution	Large 2.4" colour LCD screen
Coloured screen	No	Yes
Foetal heart rate display (FHR)	Yes	Yes
Foetal heart rate waveform display	No	Yes
Pulse indicator	Yes	Yes
Probe	2MHz	2MHz
Measuring range	50bpm - 240bpm	50bpm - 240bpm
Accuracy	+/- 2%	+/- 2%
Dimensions (L x W x D)	128mm x 75mm x 26mm	132mm x 66mm x 27mm
Weight (without batteries)	190g (without batteries)	200g (without batteries)
Waterproof	IPX7 (probe only) & IPX2	IPX7 (probe only)
Auto power off	Yes	Yes
Built in speaker	Yes	Yes
Software compatible	Yes	Yes
Standards	IEC60601-1, 1 6C60601-1-2, IEC60601-2-2	EN60601-1, EN60601-1-2
Battery life	5 hours continuous use	4 hours continuous use
Low battery indicator	Yes	Yes
Supplied with batteries	Yes (2 x 1.5V AA)	Yes (2 x 1.5V AA)
Supplied with gel	Yes	Yes
Supplied with carry case	Yes	Yes
<b>NHS SC Band 1 Price (1-9 units)</b>	<b>£79.95</b>	<b>£144.95</b>
<b>NHS SC Band 2 Price (10+ units)</b>	<b>£71.95</b>	<b>£130.50</b>

Timesco are the UK's premier manufacturer of medical devices for Anaesthesia, Surgery, Podiatry and Primary Care since 1964.

Timesco Healthcare Ltd Timesco House 3 Carnival Park  
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 710

E [info@timesco.com](mailto:info@timesco.com) | [timesco.com](http://timesco.com)

Ref: B.GM.014/2 E & O.E. Timesco Healthcare Ltd. 09/22. Prices include VAT and correct at time of issue.

**TIMESCO**  
INSTRUMENTS FOR LIFE