



CALLISTO LARYNGOSCOPE RANGE

INTRODUCING V2 & MINI HANDLES

READY TO USE • FIBRE OPTIC • SINGLE USE • EXCEPTIONAL ILLUMINATION • LED • SETS

TIMESCO
INSTRUMENTS FOR LIFE

Callisto is the U.K's No.1 single use laryngoscope system. It has been designed to offer a true low cost alternative to reusable blades and handles with exceptional specifications to reduce cross infection, cost, maintenance and processing costs.

PRODUCT FEATURES

METAL SPATULA

Rigid metal spatula looks and feels like a reusable fibre optic blade.

STRONG & SECURE

Metal hook bar offering exceptional strength for a single use handle and contains 2 metal ball bearings for secure fitting.

EXCEPTIONAL ILLUMINATION

Highly polished 4mm acrylic core for maximum illumination.

HIGH PERFORMANCE

Powerful bright LED bulb for excellent illumination with an average minimum output of 11,000 lux when measured 2cm from the bulb.

EXTENSIVE RANGE

Available in a comprehensive range of convenient ready to use sets, including Eclipse (tilting tip) for elevation of the epiglottis; to provide an improved view of the larynx without excessive force.



READY TO USE

Preloaded with alkaline button cell batteries, which are quick and easy to remove for safe battery disposal.

DISPOSABLE

Features a low volume plastic construction to reduce plastic waste and cost of disposal.

PATENT PROTECTED

China: ZL2015205363942


UK: 2534157





PRODUCT SPECIFICATION

CE CE marked

ISO 7376 ISO 7376 compatible

 Fibre-optic, green fitting

 Single Use

 Supplied with multilingual instructions [I.L.039/6](#)

 Latex free

RANGE INCLUDES



Blades




Handles




Sets

Supplied clinically clean and ready to use with a shelf life of:
Blades: 5 years
Handles: 3 years battery life

LARYNGOSCOPE BLADES



MACINTOSH




Product Code	Description	Measurements*	Quantity
DS.3940.150.05	Callisto Macintosh Size 0	80mm x 13mm	Box of 10
DS.3940.150.10	Callisto Macintosh Size 1	90mm x 18mm	Box of 10
DS.3940.150.15	Callisto Macintosh Size 2	100mm x 18mm	Box of 10
DS.3940.150.20	Callisto Macintosh Size 3	130mm x 18mm	Box of 10
DS.3940.150.23	Callisto Macintosh Size 3.5	144mm x 18mm	Box of 10
DS.3940.150.25	Callisto Macintosh Size 4	153mm x 18mm	Box of 10
DS.3940.150.30	Callisto Macintosh Size 5	174mm x 16mm	Box of 10

Macintosh design features a curved blade, are inserted into the vallecula, epiglottis fold, and moved laterally to expose the laryngeal cavity and vocal cords. These are the most popular blades worldwide.

MILLER



W x H




Product Code	Description	Measurements*	Quantity
DS.3940.185.03	Callisto Miller Size 00	65mm x 12mm	Box of 10
DS.3940.185.05	Callisto Miller Size 0	79mm x 12mm	Box of 10
DS.3940.185.10	Callisto Miller Size 1	105mm x 13mm	Box of 10
DS.3940.185.15	Callisto Miller Size 2	158mm x 15mm	Box of 10
DS.3940.185.20	Callisto Miller Size 3	197mm x 16mm	Box of 10
DS.3940.185.25	Callisto Miller Size 4	208mm x 18mm	Box of 10

Miller designs feature a straight blade, are inserted over the epiglottis and moved laterally to expose the laryngeal cavity and vocal cords. The smaller sizes 00, 0 and 1 are preferred by clinicians for paediatric and infant use.

ROBERTSHAW



W x H




Product Code	Description	Measurements*	Quantity
DS.3940.235.05	Callisto Robertshaw Size 0	94mm x 10mm	Box of 10
DS.3940.235.10	Callisto Robertshaw Size 1	112mm x 12mm	Box of 10

Due to the low profile on the Robertshaw blades, these are preferred by clinicians for paediatric and infant use.

SEWARD



W x H



Product Code	Description	Measurements*	Quantity
DS.3940.240.10	Callisto Seward Size 1	107mm x 13mm	Box of 10

Due to the low profile on the Seward blades, these are preferred by clinicians for paediatric and infant use.

In this brochure, only our Callisto single use fibre optic light guide laryngoscope systems are illustrated. We also manufacture the reusable laryngoscope systems, [Optima CLX](#), with removable fibre optic light cores, [Sirius](#) with fixed fibre optic light cores and Orion our standard light systems. Please contact us or refer to our website & separate publications for more information on these systems.

LARYNGOSCOPE BLADES

Highly polished 4mm acrylic core for maximum illumination.

Supplied in boxes of 10 when sold individually or boxes of 5 when sold as a set.

Rigid metal spatula looks and feels like a reusable fibre optic blade.



No touch hook prevents the blade from touching the handle.

Green ISO 7376 fitting plastic block compatible with green fitting single use and reusable handles.

Single use logo clearly marked on the product.

Tilting tip for elevation of the epiglottis, providing an improved view of the larynx, without excessive force.

2 metal ball bearings for secure fitting.

CALLISTO ECLIPSE BLADES

Product Code	Description	Measurements*	Quantity
DS.3940.175.20	Callisto Eclipse Blade Size 3	130mm x 18mm	Box of 10
DS.3940.175.25	Callisto Eclipse Blade Size 4	150mm x 18mm	Box of 10

*Length x Height. Height is measured from the top of the block to the top of the blade. Measurements rounded to the nearest millimetre.

Introducing 2 new sizes to the Callisto preloaded single use laryngoscope handle range. A true low cost alternative to reusable systems, the handles reduce the risk of cross infection & eliminate the cost of reprocessing. They provide exceptional strength, bright illumination from a powerful LED bulb and are made from a low volume plastic construction, producing less waste for a greener system.

Metal hook bar offering exceptional strength for a single use handle.

All the handles are ISO 7376 fitting and compatible with green fitting blades.

New lower volume plastic construction reducing plastic waste and cost of disposal.

Preloaded with alkaline button cell batteries, which are quick and easy to remove for safe battery disposal.

Available individually or with compatible blades as convenient, ready to use sets.

Powerful, bright LED bulb for excellent illumination with an average minimum output of 36,000 lux for the V2 handle or 11,000 lux for the Mini handle when measured 2cm from the bulb.

We have introduced 2 new sizes to the Callisto handle range, as we understand that one size does not fit all.

Clearly marked as single use on the base of the handle.



All handles are patent protected
China: ZL2015205363942
UK: 2534157

New sizes available



Adult V2

Mini

Product Code	Description	Measurements*	Quantity
DS.3940.450.10	Callisto Adult V2 Preloaded Handle	120mm x 39mm	Box of 25
DS.3940.490.10	Callisto Mini (Paediatric) Preloaded Handle	120mm x 26mm	Box of 10

*Height x diameter. Measurements rounded to the nearest millimetre.

CALLISTO BLADE & ADULT V2 LED PRELOADED HANDLE SETS



CALLISTO MACINTOSH BLADE & ADULT V2 LED PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.950.05	Callisto Blade Macintosh 0 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.10	Callisto Blade Macintosh 1 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.15	Callisto Blade Macintosh 2 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.20	Callisto Blade Macintosh 3 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.23	Callisto Blade Macintosh 3.5 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.25	Callisto Blade Macintosh 4 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.27	Callisto Blade Macintosh 5 & Adult V2 LED Preloaded Handle	Box of 10



CALLISTO MILLER BLADE & ADULT V2 LED PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.950.28	Callisto Blade Miller 00 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.30	Callisto Blade Miller 0 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.35	Callisto Blade Miller 1 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.40	Callisto Blade Miller 2 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.45	Callisto Blade Miller 3 & Adult V2 LED Preloaded Handle	Box of 10



CALLISTO ECLIPSE (TILTING TIP) BLADE & ADULT V2 LED PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.953.20	Callisto Blade Eclipse Macintosh 3 & Adult V2 LED Preloaded Handle	Box of 5
DS.3950.953.25	Callisto Blade Eclipse Macintosh 4 & Adult V2 LED Preloaded Handle	Box of 5

CALLISTO BLADE & MINI (PAEDIATRIC) PRELOADED HANDLE SETS



CALLISTO MACINTOSH BLADE & MINI (PAEDIATRIC) PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.955.05	Callisto Blade Macintosh 0 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.10	Callisto Blade Macintosh 1 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.15	Callisto Blade Macintosh 2 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.20	Callisto Blade Macintosh 3 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.23	Callisto Blade Macintosh 3.5 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.25	Callisto Blade Macintosh 4 & Mini (Paediatric) Preloaded Handle	Box of 10



CALLISTO MILLER BLADE & MINI (PAEDIATRIC) PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.955.28	Callisto Blade Miller 00 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.30	Callisto Blade Miller 0 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.35	Callisto Blade Miller 1 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.40	Callisto Blade Miller 2 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.45	Callisto Blade Miller 3 & Mini (Paediatric) Preloaded Handle	Box of 10

BENEFITS OF SINGLE USE

The Timesco Callisto single use sets have exceptional specifications to reduce cross infection at a cost low enough to justify single use. The sets contain single use blades and handles (preloaded with batteries) providing a ready to use disposable system to increase time efficiency and convenience.



Please refer to our website for our frequently asked questions and answers on the family of Callisto products. Visit: <https://timesco.com/single-use/>

Product Code	Description	Dimensions	Quantity
--------------	-------------	------------	----------

CALLISTO MACINTOSH LARYNGOSCOPE BLADES

DS.3940.150.05	Callisto Macintosh Size 0	80mm x 13mm	Box of 10
DS.3940.150.10	Callisto Macintosh Size 1	90mm x 18mm	Box of 10
DS.3940.150.15	Callisto Macintosh Size 2	100mm x 18mm	Box of 10
DS.3940.150.20	Callisto Macintosh Size 3	130mm x 18mm	Box of 10
DS.3940.150.23	Callisto Macintosh Size 3.5	144mm x 18mm	Box of 10
DS.3940.150.25	Callisto Macintosh Size 4	153mm x 18mm	Box of 10
DS.3940.150.30	Callisto Macintosh Size 5	174mm x 16mm	Box of 10

CALLISTO MILLER LARYNGOSCOPE BLADES

DS.3940.185.03	Callisto Miller Size 00	65mm x 12mm	Box of 10
DS.3940.185.05	Callisto Miller Size 0	79mm x 12mm	Box of 10
DS.3940.185.10	Callisto Miller Size 1	105mm x 13mm	Box of 10
DS.3940.185.15	Callisto Miller Size 2	158mm x 15mm	Box of 10
DS.3940.185.20	Callisto Miller Size 3	197mm x 16mm	Box of 10
DS.3940.185.25	Callisto Miller Size 4	208mm x 18mm	Box of 10

CALLISTO ROBERTSHAW LARYNGOSCOPE BLADES

DS.3940.235.05	Callisto Robertshaw Size 0	94mm x 10mm	Box of 10
DS.3940.235.10	Callisto Robertshaw Size 1	112mm x 12mm	Box of 10

CALLISTO SEWARD LARYNGOSCOPE BLADES

DS.3940.240.10	Callisto Seward Size 1	107mm x 13mm	Box of 10
----------------	------------------------	--------------	-----------

CALLISTO ECLIPSE LARYNGOSCOPE BLADES

DS.3940.175.20	Callisto Eclipse Blade Size 3	130mm x 18mm	Box of 10
DS.3940.175.25	Callisto Eclipse Blade Size 4	150mm x 18mm	Box of 10

CALLISTO LARYNGOSCOPE HANDLES

DS.3940.450.10	Callisto Adult V2 Preloaded Handle	120mm x 39mm	Box of 25
DS.3940.490.10	Callisto Mini (Paediatric) Preloaded Handle	120mm x 26mm	Box of 10

Length x Height. Height is measured from the top of the block to the top of the blade. Measurements rounded to the nearest millimetre.

Product Code	Description	Quantity
CALLISTO MACINTOSH BLADE & ADULT V2 LED PRELOADED HANDLE SET		
DS.3950.950.05	Callisto Blade Macintosh 0 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.10	Callisto Blade Macintosh 1 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.15	Callisto Blade Macintosh 2 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.20	Callisto Blade Macintosh 3 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.23	Callisto Blade Macintosh 3.5 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.25	Callisto Blade Macintosh 4 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.27	Callisto Blade Macintosh 5 & Adult V2 LED Preloaded Handle	Box of 10
CALLISTO MILLER BLADE & ADULT V2 LED PRELOADED HANDLE SET		
DS.3950.950.28	Callisto Blade Miller 00 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.30	Callisto Blade Miller 0 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.35	Callisto Blade Miller 1 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.40	Callisto Blade Miller 2 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.45	Callisto Blade Miller 3 & Adult V2 LED Preloaded Handle	Box of 10
CALLISTO ECLIPSE (TILTING TIP) BLADE & ADULT V2 LED PRELOADED HANDLE SET		
DS.3950.953.20	Callisto Blade Eclipse Macintosh 3 & Adult V2 LED Preloaded Handle	Box of 5
DS.3950.953.25	Callisto Blade Eclipse Macintosh 4 & Adult V2 LED Preloaded Handle	Box of 5
CALLISTO MACINTOSH BLADE & MINI (PAEDIATRIC) PRELOADED HANDLE SET		
DS.3950.955.05	Callisto Blade Macintosh 0 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.10	Callisto Blade Macintosh 1 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.15	Callisto Blade Macintosh 2 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.20	Callisto Blade Macintosh 3 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.23	Callisto Blade Macintosh 3.5 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.25	Callisto Blade Macintosh 4 & Mini (Paediatric) Preloaded Handle	Box of 10
CALLISTO MILLER BLADE & MINI (PAEDIATRIC) PRELOADED HANDLE SET		
DS.3950.955.28	Callisto Blade Miller 00 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.30	Callisto Blade Miller 0 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.35	Callisto Blade Miller 1 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.40	Callisto Blade Miller 2 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.45	Callisto Blade Miller 3 & Mini (Paediatric) Preloaded Handle	Box of 10

Contact Us

find out more today

For over 55 years' Timesco have been at the forefront of Medical product design and innovation. As a manufacturer and supplier of laryngoscopes and diagnostic products, we have an unrivalled portfolio supplying both the UK directly and internationally through distributors.

As a family owned business, with headquarters in Basildon England, we strive to deliver products made to an exceptionally high standard; exceeding internationally recognised quality standards and customer expectations. We are recognised as one the most progressive and innovative companies within the industry and have a dedicated product development team which focuses on the changing needs of our clients and advances in clinical practice.

Contact Timesco today for further information or to request product samples.

t +44 (0)1268 297 700

e sales@timesco.com

W timesco.com



ISO 13485



FDA Manufacturers
License (USA)
(Registration Number: 8021788)



NMPA
CHINA



Ministry of
Health

Classification: TIMESCO products are classified under the Harmonised Tariff: 90-1890-8400.

Timesco Healthcare Ltd, Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

TIMESCO
INSTRUMENTS FOR LIFE

T +44 (0) 1268 297 700 | **E** info@timesco.com | timesco.com

Ref: B.L.044/4 E & O.E. Timesco Healthcare Ltd. 5/22

