



CALLISTO LARYNGOSCOPE RANGE

INTRODUCING V2 & MINI HANDLES

READY TO USE • FIBRE OPTIC • SINGLE-USE • EXCEPTIONAL ILLUMINATION • LED • SETS

TIMESCO

Callisto is the U.K's No.1 single-use laryngoscope system. It has been designed to offer a true low cost alternative to reusable blades and handles with exceptional specifications to reduce cross infection, cost, maintenance and processing costs.

10 IMPORTANT REASONS TO USE SINGLE-USE

REDUCE CROSS INFECTION

Single-use products can reduce the risk of cross contamination to the patient, over reusable products that may not be effectively cleaned.



COST SAVINGS

Eliminates the cost of reprocessing products.



BRIGHT LIGHT, EVERY TIME

On average, 3x brighter than recommended in the ISO 7376 standard.



STRENGTH

The Callisto single-use laryngoscope system exceeds the 15kg strength test recommended in the ISO 7376 standard.



RECYCLING OPTIONS

The range can be recycled to help reduce the impact on the environment.



READY TO USE

Simply open, use & dispose.



LOW CHEMICAL USE

There are fewer chemicals used due to eliminating the reprocessing process.



NO WAITING

No need to wait for products to be received after the cleaning process.



SAVE TIME

No time spent transporting, cleaning & reassembling products.



PRELOADED BATTERIES


Timesco's Callisto handles contain new small button cell batteries, with a 5 year shelflife, that are easy to dispose of and mercury free.




PRODUCT SPECIFICATION

 CE marked

ISO 7376 ISO 7376 compatible

 Fibre-optic, green fitting

 Single-use

 Supplied with multilingual instructions I.L.039

 Latex free

RANGE INCLUDES



Blades



Handles



Sets

Supplied clinically clean and ready to use with a shelf life of:
Blades: 5 years
Handles: 5 years battery life

LARYNGOSCOPE BLADES

Rounded blade tip replicates the reusable Sirius fibre optic blade.

Highly polished acrylic core for maximum illumination.

Rigid, low profile, metal spatula looks and replicates a reusable Sirius fibre optic blade.

Single-use logo clearly marked on the back of the blade.

Green ISO 7376 fitting plastic block compatible with green fitting single-use and reusable handles.

'No Touch' block prevents the tip of the blade coming into contact with the handle


Supplied in boxes of 10.

Clearly stamped LOT number for tracability.

2 metal ball bearings for secure fitting.

Product description clearly marked on the back of the blade.


MACINTOSH



Product Code	Description	Measurements*	Quantity
DS.3940.150.00	Callisto Macintosh Size 00	65mm x 6.9mm	Box of 10
DS.3940.150.05	Callisto Macintosh Size 0	80mm x 6.8mm	Box of 10
DS.3940.150.10	Callisto Macintosh Size 1	9.5mm x 6.8mm	Box of 10
DS.3940.150.15	Callisto Macintosh Size 2	115mm x 6.8mm	Box of 10
DS.3940.150.20	Callisto Macintosh Size 3	134mm x 6.8mm	Box of 10
DS.3940.150.23	Callisto Macintosh Size 3.5	145mm x 6.8mm	Box of 10
DS.3940.150.25	Callisto Macintosh Size 4	155mm x 6.8mm	Box of 10
DS.3940.150.30	Callisto Macintosh Size 5	175mm x 6.8mm	Box of 10

Macintosh design features a curved blade, are inserted into the vallecula, epiglottis fold, and moved laterally to expose the laryngeal cavity and vocal cords. These are the most popular blades worldwide.

MILLER




LxH

Product Code	Description	Measurements*	Quantity
DS.3940.185.02	Callisto Miller Size 000	60mm x 12mm	Box of 10
DS.3940.185.03	Callisto Miller Size 00	65mm x 12mm	Box of 10
DS.3940.185.05	Callisto Miller Size 0	78mm x 12mm	Box of 10
DS.3940.185.10	Callisto Miller Size 1	105mm x 12mm	Box of 10
DS.3940.185.13	Callisto Miller Size 1.5	105mm x 12mm	Box of 10
DS.3940.185.15	Callisto Miller Size 2	157mm x 13.5mm	Box of 10
DS.3940.185.20	Callisto Miller Size 3	196mm x 13.5mm	Box of 10
DS.3940.185.25	Callisto Miller Size 4	208mm x 15.5mm	Box of 10

Miller designs feature a straight blade, are inserted over the epiglottis and moved laterally to expose the laryngeal cavity and vocal cords. The smaller sizes 000, 00 and 0 are preferred by clinicians for paediatric and infant use.

ROBERTSHAW




LxH

Product Code	Description	Measurements*	Quantity
DS.3940.235.05	Callisto Robertshaw Size 0	94mm x 11.7mm	Box of 10
DS.3940.235.10	Callisto Robertshaw Size 1	112mm x 13.7mm	Box of 10

Due to the low profile on the Robertshaw blades, these are preferred by clinicians for paediatric and infant use.

SEWARD



LxH

Product Code	Description	Measurements*	Quantity
DS.3940.240.10	Callisto Seward Size 1	107mm x 10.3mm	Box of 10

Due to the low profile on the Seward blades, these are preferred by clinicians for paediatric and infant use.

In this brochure, only our Callisto single-use fibre optic light guide laryngoscope systems are illustrated. We also manufacture the reusable laryngoscope systems, [Optima CLX](#), with removable fibre optic light cores, [Sirius](#) with fixed fibre optic light cores and Orion our standard light systems. Please contact us or refer to our website & separate publications for more information on these systems.

LARYNGOSCOPE BLADES

Highly polished acrylic core for maximum illumination.

Supplied in boxes of 10 when sold individually or boxes of 5 when sold as a set.

Rigid metal spatula looks and feels like a reusable fibre optic blade.

Single-use logo clearly marked on the product.



Tilting tip for elevation of the epiglottis, providing an improved view of the larynx, without excessive force.

Green ISO 7376 fitting plastic block compatible with green fitting single-use and reusable handles.

2 metal ball bearings for secure fitting.



CALLISTO ECLIPSE BLADES

Product Code	Description	Measurements*	Quantity
DS.3940.175.20	Callisto Eclipse Blade Size 3	134mm x 14mm	Box of 10
DS.3940.175.25	Callisto Eclipse Blade Size 4	155mm x 14mm	Box of 10

*Length x Height. Height is measured from the top of the block to the top of the blade. Measurements rounded to the nearest millimetre.

Introducing 2 new sizes to the Callisto preloaded single-use laryngoscope handle range. A true low cost alternative to reusable systems, the handles reduce the risk of cross infection & eliminate the cost of reprocessing. They provide exceptional strength, bright illumination from a powerful LED bulb and are made from a low volume plastic construction, producing less waste for a greener system.



Metal hook bar offering exceptional strength for a single-use handle.

All the handles are ISO 7376 fitting and compatible with green fitting blades.

New lower volume plastic clam shell construction reduces plastic waste and cost of disposal, whilst increasing handle strength.

Preloaded with alkaline button cell batteries that have a 5 year shelflife. Easily dispose of the batteries in 3 quick and easy steps.
1. Twist. 2. Drop. 3. Dispose.



New sizes available



Powerful, bright LED bulb for excellent illumination with an average minimum output of 36,000 lux for the V2 handle or 11,000 lux for the Mini handle when measured 2cm from the bulb.

We've expanded the Callisto handle range to encompass 3 different sizes, providing clinicians with a selection of handles to meet their specific requirements.

Clearly marked with product code, CE mark and single-use symbol on the base.



All handles are patent protected
China: ZL2015205363942
UK: 2534157

Available individually or with compatible blades as convenient, ready to use sets.



New

Adult V2

New

Mini

Compact

CALLISTO PRELOADED HANDLES

Product Code	Description	Measurements*	Quantity
DS.3940.460.10	Callisto Compact Preloaded Handle	120mm x 29mm	Box of 25
DS.3940.450.10	Callisto Adult V2 Preloaded Handle	119mm x 33mm	Box of 25
DS.3940.490.10	Callisto Mini (Paediatric) Preloaded Handle	120mm x 22mm	Box of 10

*Height x diameter. Measurements rounded to the nearest millimetre.

CALLISTO BLADE & COMPACT PRELOADED HANDLE SETS



CALLISTO MACINTOSH BLADE & COMPACT PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.952.00	Callisto Blade Macintosh 00 & Compact Preloaded Handle	Box of 10
DS.3950.952.05	Callisto Blade Macintosh 0 & Compact Preloaded Handle	Box of 10
DS.3950.952.10	Callisto Blade Macintosh 1 & Compact Preloaded Handle	Box of 10
DS.3950.952.15	Callisto Blade Macintosh 2 & Compact Preloaded Handle	Box of 10
DS.3950.952.20	Callisto Blade Macintosh 3 & Compact Preloaded Handle	Box of 10
DS.3950.952.25	Callisto Blade Macintosh 4 & Compact Preloaded Handle	Box of 10
DS.3950.952.65	Callisto Blade Macintosh 5 & Compact Preloaded Handle	Box of 10



CALLISTO MILLER BLADE & COMPACT PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.952.27	Callisto Blade Miller 000 & Compact Preloaded Handle	Box of 10
DS.3950.952.28	Callisto Blade Miller 00 & Compact Preloaded Handle	Box of 10
DS.3950.952.30	Callisto Blade Miller 0 & Compact Preloaded Handle	Box of 10
DS.3950.952.35	Callisto Blade Miller 1 & Compact Preloaded Handle	Box of 10
DS.3950.952.40	Callisto Blade Miller 2 & Compact Preloaded Handle	Box of 10
DS.3950.952.45	Callisto Blade Miller 3 & Compact Preloaded Handle	Box of 10



CALLISTO ROBERTSHAW BLADE & COMPACT PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.952.55	Callisto Blade Robertshaw 0 & Compact Preloaded Handle	Box of 10
DS.3950.952.60	Callisto Blade Robertshaw 1 & Compact Preloaded Handle	Box of 10

CALLISTO BLADE & ADULT V2 LED PRELOADED HANDLE SETS



CALLISTO MACINTOSH BLADE & ADULT V2 LED PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.950.00	Callisto Blade Macintosh 00 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.05	Callisto Blade Macintosh 0 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.10	Callisto Blade Macintosh 1 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.15	Callisto Blade Macintosh 2 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.20	Callisto Blade Macintosh 3 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.23	Callisto Blade Macintosh 3.5 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.25	Callisto Blade Macintosh 4 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.27	Callisto Blade Macintosh 5 & Adult V2 LED Preloaded Handle	Box of 10



CALLISTO MILLER BLADE & ADULT V2 LED PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.950.65	Callisto Blade Miller 000 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.28	Callisto Blade Miller 00 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.30	Callisto Blade Miller 0 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.35	Callisto Blade Miller 1 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.40	Callisto Blade Miller 2 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.45	Callisto Blade Miller 3 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.50	Callisto Blade Miller 4 & Adult V2 LED Preloaded Handle	Box of 10



CALLISTO ECLIPSE (TILTING TIP) BLADE & ADULT V2 LED PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.953.20	Callisto Blade Eclipse Macintosh 3 & Adult V2 LED Preloaded Handle	Box of 5
DS.3950.953.25	Callisto Blade Eclipse Macintosh 4 & Adult V2 LED Preloaded Handle	Box of 5

CALLISTO BLADE & MINI (PAEDIATRIC) PRELOADED HANDLE SETS



CALLISTO MACINTOSH BLADE & MINI (PAEDIATRIC) PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.955.00	Callisto Blade Macintosh 00 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.05	Callisto Blade Macintosh 0 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.10	Callisto Blade Macintosh 1 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.15	Callisto Blade Macintosh 2 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.20	Callisto Blade Macintosh 3 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.23	Callisto Blade Macintosh 3.5 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.25	Callisto Blade Macintosh 4 & Mini (Paediatric) Preloaded Handle	Box of 10



CALLISTO MILLER BLADE & MINI (PAEDIATRIC) PRELOADED HANDLE SET

Product Code	Description	Quantity
DS.3950.955.27	Callisto Blade Miller 000 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.28	Callisto Blade Miller 00 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.30	Callisto Blade Miller 0 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.35	Callisto Blade Miller 1 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.40	Callisto Blade Miller 2 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.45	Callisto Blade Miller 3 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.50	Callisto Blade Miller 4 & Mini (Paediatric) Preloaded Handle	Box of 10

BENEFITS OF SINGLE-USE

The Timesco Callisto single-use sets have exceptional specifications to reduce cross infection at a cost low enough to justify single-use. The sets contain single-use blades and handles (preloaded with batteries) providing a ready to use disposable system to increase time efficiency and convenience.



Please refer to our website for our frequently asked questions and answers on the family of Callisto products. Visit: <https://timesco.com/single-use/>

Product Code	Description	Dimensions	Quantity
CALLISTO MACINTOSH LARYNGOSCOPE BLADES			
DS.3940.150.00	Callisto Macintosh Size 00	65mm x 6.9mm	Box of 10
DS.3940.150.05	Callisto Macintosh Size 0	80mm x 6.8mm	Box of 10
DS.3940.150.10	Callisto Macintosh Size 1	9.5mm x 6.8mm	Box of 10
DS.3940.150.15	Callisto Macintosh Size 2	115mm x 6.8mm	Box of 10
DS.3940.150.20	Callisto Macintosh Size 3	134mm x 6.8mm	Box of 10
DS.3940.150.23	Callisto Macintosh Size 3.5	145mm x 6.8mm	Box of 10
DS.3940.150.25	Callisto Macintosh Size 4	155mm x 6.8mm	Box of 10
DS.3940.150.30	Callisto Macintosh Size 5	175mm x 6.8mm	Box of 10
CALLISTO MILLER LARYNGOSCOPE BLADES			
DS.3940.185.02	Callisto Miller Size 000	60mm x 12mm	Box of 10
DS.3940.185.03	Callisto Miller Size 00	65mm x 12mm	Box of 10
DS.3940.185.05	Callisto Miller Size 0	78mm x 12mm	Box of 10
DS.3940.185.10	Callisto Miller Size 1	105mm x 12mm	Box of 10
DS.3940.185.13	Callisto Miller Size 1.5	132mm x 13.5mm	Box of 10
DS.3940.185.15	Callisto Miller Size 2	157mm x 13.5mm	Box of 10
DS.3940.185.20	Callisto Miller Size 3	196mm x 13.5mm	Box of 10
DS.3940.185.25	Callisto Miller Size 4	208mm x 15.5mm	Box of 10
CALLISTO ROBERTSHAW LARYNGOSCOPE BLADES			
DS.3940.235.05	Callisto Robertshaw Size 0	94mm x 11.7mm	Box of 10
DS.3940.235.10	Callisto Robertshaw Size 1	112mm x 13.7mm	Box of 10
CALLISTO SEWARD LARYNGOSCOPE BLADES			
DS.3940.240.10	Callisto Seward Size 1	107mm x 10.3mm	Box of 10
CALLISTO ECLIPSE LARYNGOSCOPE BLADES			
DS.3940.175.20	Callisto Eclipse Blade Size 3	134mm x 14mm	Box of 10
DS.3940.175.25	Callisto Eclipse Blade Size 4	155mm x 14mm	Box of 10
CALLISTO LARYNGOSCOPE HANDLES			
DS.3940.460.10	Callisto Compact Preloaded Handle	120mm x 29mm	Box of 25
DS.3940.450.10	Callisto Adult V2 Preloaded Handle	119mm x 33mm	Box of 25
DS.3940.490.10	Callisto Mini (Paediatric) Preloaded Handle	120mm x 22mm	Box of 10

Length x Height. Height is measured from the top of the block to the top of the blade. Measurements rounded to the nearest millimetre.

PRODUCT LIST (CONTINUED)

Product Code	Description	Quantity
CALLISTO MACINTOSH BLADE & COMPACT PRELOADED HANDLE SET		
DS.3950.952.00	Callisto Blade Macintosh 00 & Compact Preloaded Handle	Box of 10
DS.3950.952.05	Callisto Blade Macintosh 0 & Compact Preloaded Handle	Box of 10
DS.3950.952.10	Callisto Blade Macintosh 1 & Compact Preloaded Handle	Box of 10
DS.3950.952.15	Callisto Blade Macintosh 2 & Compact Preloaded Handle	Box of 10
DS.3950.952.20	Callisto Blade Macintosh 3 & Compact Preloaded Handle	Box of 10
DS.3950.952.25	Callisto Blade Macintosh 4 & Compact Preloaded Handle	Box of 10
DS.3950.952.65	Callisto Blade Macintosh 5 & Compact Preloaded Handle	Box of 10
CALLISTO MILLER BLADE & COMPACT PRELOADED HANDLE SET		
DS.3950.952.27	Callisto Blade Miller 000 & Compact Preloaded Handle	Box of 10
DS.3950.952.28	Callisto Blade Miller 00 & Compact Preloaded Handle	Box of 10
DS.3950.952.30	Callisto Blade Miller 0 & Compact Preloaded Handle	Box of 10
DS.3950.952.35	Callisto Blade Miller 1 & Compact Preloaded Handle	Box of 10
DS.3950.952.40	Callisto Blade Miller 2 & Compact Preloaded Handle	Box of 10
DS.3950.952.45	Callisto Blade Miller 3 & Compact Preloaded Handle	Box of 10
CALLISTO ROBERTSHAW BLADE & COMPACT PRELOADED HANDLE SET		
DS.3950.952.55	Callisto Blade Robershaw 0 & Compact Preloaded Handle	Box of 10
DS.3950.952.60	Callisto Blade Robershaw 1 & Compact Preloaded Handle	Box of 10
CALLISTO MACINTOSH BLADE & ADULT V2 LED PRELOADED HANDLE SET		
DS.3950.950.00	Callisto Blade Macintosh 00 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.05	Callisto Blade Macintosh 0 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.10	Callisto Blade Macintosh 1 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.15	Callisto Blade Macintosh 2 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.20	Callisto Blade Macintosh 3 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.23	Callisto Blade Macintosh 3.5 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.25	Callisto Blade Macintosh 4 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.27	Callisto Blade Macintosh 5 & Adult V2 LED Preloaded Handle	Box of 10
CALLISTO MILLER BLADE & ADULT V2 LED PRELOADED HANDLE SET		
DS.3950.950.65	Callisto Blade Miller 000 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.28	Callisto Blade Miller 00 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.30	Callisto Blade Miller 0 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.35	Callisto Blade Miller 1 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.40	Callisto Blade Miller 2 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.45	Callisto Blade Miller 3 & Adult V2 LED Preloaded Handle	Box of 10
DS.3950.950.50	Callisto Blade Miller 4 & Adult V2 LED Preloaded Handle	Box of 10

Product Code	Description	Quantity
--------------	-------------	----------

CALLISTO ECLIPSE (TILTING TIP) BLADE & ADULT V2 LED PRELOADED HANDLE SET

DS.3950.953.20	Callisto Blade Eclipse Macintosh 3 & Adult V2 LED Preloaded Handle	Box of 5
DS.3950.953.25	Callisto Blade Eclipse Macintosh 4 & Adult V2 LED Preloaded Handle	Box of 5

CALLISTO MACINTOSH BLADE & MINI (PAEDIATRIC) PRELOADED HANDLE SET

DS.3950.955.00	Callisto Blade Macintosh 00 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.05	Callisto Blade Macintosh 0 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.10	Callisto Blade Macintosh 1 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.15	Callisto Blade Macintosh 2 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.20	Callisto Blade Macintosh 3 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.23	Callisto Blade Macintosh 3.5 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.25	Callisto Blade Macintosh 4 & Mini (Paediatric) Preloaded Handle	Box of 10

CALLISTO MILLER BLADE & MINI (PAEDIATRIC) PRELOADED HANDLE SET

DS.3950.955.27	Callisto Blade Miller 000 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.28	Callisto Blade Miller 00 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.30	Callisto Blade Miller 0 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.35	Callisto Blade Miller 1 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.40	Callisto Blade Miller 2 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.45	Callisto Blade Miller 3 & Mini (Paediatric) Preloaded Handle	Box of 10
DS.3950.955.50	Callisto Blade Miller 4 & Mini (Paediatric) Preloaded Handle	Box of 10



ISO 13485



FDA Manufacturers
License (USA)
(Registration Number: 8021788)



NMPA
CHINA



Ministry of Health
Ministry of Health

Classification: TIMESCO products are classified under the Harmonised Tariff: 90-1890-8400.

Timesco Healthcare Ltd., 3 Carnival Park, Carnival Close,
Basildon, Essex, SS14 3WN, United Kingdom
T: +44 (0)1268 297 700 E: info@timesco.com W: timesco.com



Ref: B.L.044/9 FEB24 E & O.E. Timesco Healthcare Ltd. May 2023