

callisto
TIMESCO™

The Callisto preloaded single use laryngoscope handle provides bright illumination from a powerful LED bulb for optimal visualisation.

Offering a true low cost alternative to reusable systems; the handles provide exceptional strength whilst reducing the risk of cross infection & eliminating the cost of reprocessing. Available in both medium and compact size.



CE marked

**ISO
7376**

ISO 7376
compatible



Single Use



Supplied with
multilingual
instructions



Latex Free



Compatible
blades are
colour coded

SPECIFICATION

FEATURES

Metal hook bar offering exceptional strength for a single use handle.

This handle is ISO 7376 fitting and compatible with green fitting blades.

Preloaded with alkaline button cell batteries, which are quick and easy to remove for safe battery disposal.

Powerful, bright LED bulb for excellent illumination with a average minimum output of 11,000 lux when measured 2cm from the bulb.

Supplied in boxes of 25.

Low volume plastic construction reducing plastic waste and cost of disposal.



Clearly marked as single use on the base of the handle.



Medium



Compact

Product Code	Description	Measurements*	NPC Code
DS.3940.480.10	Callisto Preloaded Handle, Medium	139mm x 33mm	FSM1747
DS.3940.460.10	Callisto Preloaded Handle, Compact	120mm x 29mm	FSM3228

Timesco Healthcare Ltd Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800

E info@timesco.com | [timesco.com](https://www.timesco.com)

L.L.055/2 E & O.E Timesco Nov 2019.

*Height x Diameter. Measurements rounded to the nearest millimeter.