

# easy bright

TIMESCO™

The Timesco Easy Bright conventional single use sets offer a convenient, compact, lightweight system. They are supplied ready to use with a preloaded handle and rigid metal blade.



CE marked

**ISO  
7376**

ISO 7376  
compatible



Single Use



Supplied with  
multilingual  
instructions



Latex Free



Conventional,  
standard light  
system

## SPECIFICATION

## FEATURES

Bright cool touch LED bulb.

Rigid metal spatula looks and feels like a reusable blade.

Standard light ISO 7376  
handle and blade.

Conveniently supplied connected  
as a ready to use set for testing in  
the pouch.

Preloaded with button cells  
for convenience which are  
removable for disposal if  
required.

Supplied in packs of 10  
in easy to open pouches.

A lightweight, compact set  
weighing only 93 grams\*.

(\*Macintosh 3 set)





Product Code	Description	Dimensions of Handle <sup>1</sup>	Dimensions of Blade <sup>2</sup>
DS.2960.952.05	Easy Bright single use blade and preloaded handle set - Macintosh 0	118mm x 25mm	85mm x 14mm
DS.2960.952.10	Easy Bright single use blade and preloaded handle set - Macintosh 1	118mm x 25mm	90mm x 20mm
DS.2960.952.15	Easy Bright single use blade and preloaded handle set - Macintosh 2	118mm x 25mm	100mm x 20mm
DS.2960.952.20	Easy Bright single use blade and preloaded handle set - Macintosh 3	118mm x 25mm	130mm x 20mm
DS.2960.952.25	Easy Bright single use blade and preloaded handle set - Macintosh 4	118mm x 25mm	153mm x 20mm
DS.2960.952.28	Easy Bright single use blade and preloaded handle set - Miller 00	118mm x 25mm	65mm x 14mm
DS.2960.952.30	Easy Bright single use blade and preloaded handle set - Miller 0	118mm x 25mm	79mm x 14mm
DS.2960.952.35	Easy Bright single use blade and preloaded handle set - Miller 1	118mm x 25mm	105mm x 15mm
DS.2960.952.40	Easy Bright single use blade and preloaded handle set - Miller 2	118mm x 25mm	153mm x 17mm
DS.2960.952.45	Easy Bright single use blade and preloaded handle set - Miller 3	118mm x 25mm	197mm x 17mm
DS.2960.952.50	Easy Bright single use blade and preloaded handle set - Miller 4	118mm x 25mm	208mm x 20mm

Timesco Healthcare Ltd Timesco House 3 Carnival Park  
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800

E info@timesco.com | [timesco.com](http://timesco.com)