

easy bright ECLIPSE

TIMESCO™

The conventional Timesco Easy Bright Eclipse sets offer a single use, cost effective and ready to use solution to assist with difficult intubations.

The sets comprise of a lightweight preloaded handle and rigid metal Eclipse blade.



CE marked

**ISO
7376**

ISO 7376
compatible



Single Use



Supplied with
multilingual
instructions



Latex Free



Conventional,
standard light
system

SPECIFICATION

Tilting tip for elevation of the epiglottis, providing an improved view of the larynx without excessive force.

Bright cool touch LED bulb.

Rigid metal spatula looks and feels like a reusable blade.

Standard light ISO 7376 handle and blade.

Conveniently supplied connected as a ready to use set for testing in the pouch.

Preloaded with button cells for convenience which are removable for disposal if required.

Lever positioned at the rear of the blade to control the tilting tip.

Single use logo clearly marked on the back of the blade.

Supplied in boxes of 10.

FEATURES





Product Code	Description	Dimensions of Handle ¹	Dimensions of Blade ²
DS.2960.953.20	Easy Bright Eclipse Set - Size 3	118mm x 25mm	130 (L) x 17 mm (H)
DS.2960.953.25	Easy Bright Eclipse Set - Size 4	118mm x 25mm	153 (L) x 17 mm (H)

Timesco Healthcare Ltd Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800

E info@timesco.com | timesco.com

L.L.062/2 E & O.E Timesco Nov 2019

¹ Height x Diameter. ² Length x Height. Height is measured from the top of the block to the top of the blade. Measurements rounded to the nearest millimeter.