

# optima<sup>CLX</sup>

TIMESCO™

Optima CLX are reusable blades with a removable glass fibre light guide.

The light guide easily clicks in and out without the use of tools to allow cleaning and replacement of the core, thus reducing the cost of ownership and maintenance.



CE marked



Fully autoclaveable



Latex Free



Extensive guarantees

**ISO  
7376**

ISO 7376 compatible



Supplied with multilingual instructions



Compatible handles are colour coded

## SPECIFICATION

## FEATURES

The fibre core is positioned to offer the least obstructive view.



Clear product information etched on the blade.

3 metal ball bearings for secure fitting.

Highly polished 4mm replaceable, glass fibre core for focused illumination. It has a one piece snap in/out fibre core for easy cleaning and replacement. No tools required.



A green, fibre optic ISO 7376 fitting metal block.

Atraumatic distal end.

Fully autoclaveable metal blade with an anti-glare satin finish. Covered by a 5 year guarantee or 4,000 autoclaves, whichever is sooner. See website for full terms and conditions.



## MACINTOSH

Product Code	Description	Measurements*	NPC Code
2970.150.05	Optima CLX Macintosh Size 0	79mm x 15.5mm	FSM765
2970.150.10	Optima CLX Macintosh Size 1	92mm x 22mm	FSM766
2970.150.15	Optima CLX Macintosh Size 2	100mm x 21mm	-----
2970.150.20	Optima CLX Macintosh Size 3	131mm x 23.5mm	FSM768
2970.150.25	Optima CLX Macintosh Size 4	155mm x 25mm	FSM769



## MILLER

Product Code	Description	Measurements*	NPC Code
2970.185.03	Optima CLX Miller Size 00	66.5mm x 11.5mm	FSM773
2970.185.05	Optima CLX Miller Size 0	77mm x 11.5mm	FSM774
2970.185.10	Optima CLX Miller Size 1	101.5mm x 12mm	FSM775
2970.185.15	Optima CLX Miller Size 2	153mm x 15mm	-----
2970.185.20	Optima CLX Miller Size 3	197mm x 15mm	FSM777
2970.185.25	Optima CLX Miller Size 4	205mm x 17mm	FSM778

In addition to the Optima CLX blades above, we also supply speciality blades, including Robertshaw, Seward and Wisconsin patterns.

Spare Optima CLX cores are also available.

**Please call 01268 297 710 or visit [www.timesco.com](http://www.timesco.com) for further details.**

Timesco Healthcare Ltd Timesco House 3 Carnival Park  
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800  
E [info@timesco.com](mailto:info@timesco.com) | [timesco.com](http://timesco.com)

L.L.059/1 E & O.E Timesco Feb 2018.

\*Length x Height. Height is measured from the top of the block to the top of the blade. Measurements rounded to the nearest 0.5 millimeter.