

SINGLE USE FACE MASKS

TIMESCO™

The Timesco single use air cushion face masks reduce the risk of cross infection, provide exceptional patient comfort and are clear for optimum patient visibility. Available in a full range of sizes from neonatal to large adult.



CE Marked

ISO

ISO 10651
ISO 5356-1
ISO 10993



Single Use



Supplied with
multilingual
instructions



Latex Free

SPECIFICATION

The transparent mask allows for clear patient visibility.

Clearly marked sizing provided on each mask – available from size 0 (neonate) to size 5 (large adult).

A soft air cushion provides exceptional patient comfort.



Also available with an adjustable check valve ensuring minimum leakage and maximum patient comfort.

Single use to reduce the risk of cross infection.

Supplied with a 4 hook ring for use with a head strap.

Masks also available with a lemon scent.

FEATURES



WITH HOOK RING

Box of 50	Description	Size	NPC Code
TDM-MP-1200	Neonate	0	FDD2570
TDM-MP-1201	Infant	1	FDD2569
TDM-MP-1202	Paediatric	2	FDD2565

Box of 32	Description	Size	NPC Code
TDM-MP-1203	Small Adult	3	FDD2566
TDM-MP-1204	Medium Adult	4	FDD2567
TDM-MP-1205	Large Adult	5	FDD5153

WITH CHECK VALVE & HOOK RING

Box of 50	Description	Size	NPC Code
TDM-MP-1210	Neonate	0	FDD5154
TDM-MP-1211	Infant	1	FDD2561
TDM-MP-1212	Paediatric	2	FDD2559

Box of 32	Description	Size	NPC Code
TDM-MP-1213	Small Adult	3	FDD2564
TDM-MP-1214	Medium Adult	4	FDD2562
TDM-MP-1215	Large Adult	5	FDD5155



WITH CHECK VALVE, HOOK RING & SCENTED

Box of 50	Description	Size	NPC Code
TDM-MP-1220	Neonate	0	FDD2556
TDM-MP-1221	Infant	1	FDD5156
TDM-MP-1222	Paediatric	2	FDD5157

Box of 32	Description	Size	NPC Code
TDM-MP-1223	Small Adult	3	FDD5158
TDM-MP-1224	Medium Adult	4	FDD5159
TDM-MP-1225	Large Adult	5	FDD5160

Timesco Healthcare Ltd Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800
E info@timesco.com | timesco.com

L.L.061/2 E & O.E Timesco February 2019