

xenon TIMESCO™

The Timesco Xenon handles are available in medium, small & stubby size. They are fully autoclavable and have an extensive guarantee. They offer exceptionally bright illumination from a replaceable Xenon bulb and can be used with reusable and single use blades.



CE marked



Latex Free



Fully autoclavable

**ISO
7376**

ISO 7376 compatible



Supplied with multilingual instructions



Compatible blades are colour coded



Extensive guarantee

SPECIFICATION

A polished metal collar reflects and directs light up towards the fibre core of the blade for maximum illumination.

Clear product branding with the product code and CE located on the opposite side to the Xenon logo. The LOT number is shown on the top of the handle.

Knurled handle design provides excellent grip.

Fully autoclavable metal handle (excluding bulbs), covered by a 5 year guarantee or 4,000 autoclaves, whichever is sooner. See website for full terms and conditions.



Super-bright replaceable 2.5V Xenon bulb providing excellent, white illumination with a LUX output of 9,750. Bulbs provide up to 15 hours of use and a colour temperature of 3,050(k).

Colour coded band to clearly identify it is ISO 7376 fitting compatible with green fitting reusable and single use blades.

Powered by dry cell batteries or easily converted to a rechargeable handle by using the Timesco R Lithium Engines.

Secure screw cap for easy access to batteries.

FEATURES



Small



Medium



Stubby

Product Code	Description	Measurements*	NPC Code
3000.350.05	Xenon Handle, Small	146mm x 19mm	FSM426
3000.350.10	Xenon Handle, Medium	146mm x 29mm	FSM065
3000.358.10	Xenon Handle, Stubby	111mm x 32mm	FSM428
3000.900.20	Spare Bulbs (Pack of 6)	-----	-----

Timesco Healthcare Ltd Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

T +44 (0) 1268 297 700 | F +44 (0) 1268 297 800
E info@timesco.com | timesco.com

LL.051/2 E & O.E Timesco March 2019.

*Height x Diameter. Measurements rounded to the nearest millimeter.