



TIMESCO ^R RECHARGERS AND BATTERIES

INNOVATIVE DESIGN • LITHIUM ION • INTELLIGENT CHARGING • MULTIPLE CHARGING PORTS

TIMESCO R™

RECHARGERS

BATTERIES

An innovative, intelligent Lithium charger for use with all Timesco rechargeable engines. Both sides of the charger have end panels that can be quickly removed to expose the unique connector plates. These connector plates allow up to 3 chargers to be daisy chained together delivering power through all chargers equally only using one transformer.

We have also developed Lithium ION batteries that are fully compatible with this new charger and can be used in our existing metal laryngoscope handles; Xenon, XLED, Orion and in our Orion and Optima diagnostic handles. Batteries are available in medium C and small AA sizes.



Once engines are fully charged the recharger will switch to a low power mode to avoid using unnecessary electricity.

Only the engine is inserted into the charger, reducing the risk of cross contamination from a handle to the charger.

A fully charged Timesco R Engine Medium (C Size) will last approx. 2 hours.

The ION-R charger (single and daisy chained) takes approx. 4 hours to fully charge batteries for medium handles.

Side plates can be removed to daisy chain 3 chargers together to receive power from a single power supply.

The colour changing LED indicates the health of the battery and the charging cycle. A yellow light indicates testing and a green light indicates fully charged batteries.



PRODUCT SPECIFICATION

 CE marked

 Latex free

 Extensive guarantee



Supplied with multilingual instructions
[IL022/2](#)

89/336/EEC

Electromagnetic compatibility directive (as amended)

73/23/EEC

Low voltage directive (as amended)



TIMESCO ION-R RECHARGERS

Product Code	Description
3125.999.10	Timesco ION-R Charger (Single) Includes: Charger, power supply and both Timesco R inserts (AA and C)
3125.999.20	Timesco ION-R Charger (Dual) Includes: 2 x Chargers, 1 x power supply and 2 x each Timesco R inserts (AA and C)

Timesco ION-R Rechargers are multi voltage and are supplied with accessories for operation worldwide.



Battery for
Small Handles



Battery for
Medium Handles

Only the Timesco R Engine is inserted into the charger.

Both the small and medium Timesco R Engines can be charged in the single and double chargers.

Timesco R Engines can be fully charged in approx. 4 hours.

The colour changing LED indicates the health of the battery and the charging cycle.

Side plates can be removed to daisy chain 3 chargers together to receive power from a single power supply.

Once Timesco R Engines are fully charged the recharger will switch to a low power mode to avoid using unnecessary electricity.

TIMESCO R BATTERIES

Product Code	Description
3115.150.20	Timesco R Engine Small (AA Size) 2.5V
3115.150.50	Timesco R Engine Medium (C Size) 2.5V

Only one battery required per handle. Timesco supply small and medium rechargeable batteries to help you create a rechargeable system.

Product Code	Description	Dimensions
TIMESCO BATTERIES		
3115.150.20	Timesco R Engine Small (AA Size) 2.5V	Small (AA Size)
3115.150.50	Timesco R Engine Medium (C Size) 2.5V	Medium (C Size)
TIMESCO ION-R CHARGERS		
3125.999.10	Timesco ION-R Charger (Single) Includes: Charger, power supply and both Timesco R inserts (AA and C)	
3125.999.20	Timesco ION-R Charger (Dual) Includes: 2 x Chargers, 1 x power supply and 2 x each Timesco R inserts	

Contact Us

find out more today

For over 55 years' Timesco have been at the forefront of Medical product design and innovation. As a manufacturer and supplier of laryngoscopes and diagnostic products, we have an unrivalled portfolio supplying both the UK directly and internationally through distributors.

As a family owned business, with headquarters in Basildon England, we strive to deliver products made to an exceptionally high standard; exceeding internationally recognised quality standards and customer expectations. We are recognised as one the most progressive and innovative companies within the industry and have a dedicated product development team which focuses on the changing needs of our clients and advances in clinical practice.

Contact Timesco today for further information or to request product samples.

t +44 (0)1268 297 700

e sales@timesco.com

w timesco.com



ISO 13485



FDA Manufacturers
License (USA)
(Registration Number: 8021788)



NMPA
CHINA



Ministry of Health

Classification: TIMESCO products are classified under the Harmonised Tariff: 90-1890-8400.

Timesco Healthcare Ltd, Timesco House 3 Carnival Park
Carnival Close Basildon Essex SS14 3WN England

TIMESCO
INSTRUMENTS FOR LIFE

T +44 (0) 1268 297 700 | E info@timesco.com | timesco.com

Ref:B.L.047/4 E&O.E. Timesco Healthcare Ltd. 03/22

