



GENERAL MEDICAL

PROFESSIONAL ESSENTIALS
TIMESCO

PRODUCT RANGE	PAGE	
Wira Flexi	2	
Thermometers	3-4	
Sphygmomanometers	5-7	
Blood Pressure Monitors	8	
Stethoscopes	9-12	
Sphygmomanometer & Stethoscopes Sets & Pink Products	13	
Patella Hammers, Tuning Forks & Ring Cutters	14	
Diagnostic Sets (Otosopes/Ophthalmoscopes)	15-18	
Pen Torches	19	
NeuroLight	20	
Scissors	21	
Watches & Tongue Depressors	22	
TourniButton, TourniBand & TourniClip	23-24	
Pulse Oximeters	25-28	
Fetal Dopplers	29	
Seca Scales	30	
Nebulisers / Respiratory	31-32	
Suction Devices	33	
Resuscitation Aids	34-35	
Face Masks	36	
Airway Management, Forceps, ET Tubes & LMAs	37-38	
Laryngoscope Range	39-42	
Vide Chest Drainage Unit	43-44	

WHO WE ARE



REPUTABLE BRAND.

As a family owned business, Timesco have been at the forefront of medical product design & innovation since 1964 & are recognised as one of the most progressive & innovative companies within the industry. We are proud to manufacture & supply an unrivalled portfolio of medical products to medical professionals in Anaesthesia, Surgery, Podiatry & General Medical.



END USER SATISFACTION.

It's important that the end user is satisfied with the products that we supply. From concept through to design & manufacture, our aim is to deliver products made to exceed internationally recognised quality standards & customer expectations.



CONVENIENT LOGISTICS.


All of our products are available to order direct, alternatively we work with third party distributors including NHS Supply Chain & Bunzl. The majority of our General Medical products are available across a number of NHS Supply Chain framework agreements, & their NPC codes can be identified throughout the brochure.




SALES SUPPORT.

We have a dedicated external field sales team strategically located offering extensive coverage across the UK (Territory Manager details are available on request). We also have an internal sales customer support team based at our head office in Basildon, Essex. Our sales & support team are available Monday to Friday from 9am to 5pm.

Tel: 01268 297710 | Email: sales@timesco.com | Live chat: www.timesco.com




Breathing Tubes



Urinary Catheters



Cannulas



Chest Drains

★ NEW
PRODUCT

WIRA-FLEXI – NON-WOVEN CATHETER AND TUBE FIXATION DEVICE

STERILE | 2 |

- Strong, re-sealable hook and loop
- Hypo-allergenic & latex free
- Flexible and breathable material
- Pre-cut and ready to use

DESCRIPTION	PRODUCT CODE	NPC CODE
Wira-Flexi S (80mm x 30mm) Box of 50	C80.800	N/A
Wira-Flexi L (90mm x 35mm) Box of 50	C80.805	FEC85000
Wira-Flexi XL (150mm x 45mm) Box of 25	C80.810	FEC85001



DIGITAL DELUXE TYMPANIC EAR THERMOMETER

- Large backlit LCD: 41.0 x 20.0mm
- Measurement response time of 1 second
- Memory capacity of 10 memory readings
- °C & °F readings
- Accuracy $\pm 0.2^{\circ}\text{C}$ from 35.5°C to 42°C
- Audible beeper
- Automatic switch-off
- Low battery indication
- Battery: 2*AA
- Carry bag included
- Requires disposable probe cover (see D60.210 / FWH121)

DESCRIPTION	PRODUCT CODE	NPC CODE
Digital Delux Ear Thermometer	D60.205	FWH348
Ear Thermometer Probe Covers	D60.210	FWH121



DIGITAL ECONOMY TYMPANIC EAR THERMOMETER

- Compact & ergonomic
- Measurement response time of 1 second
- Memory capacity of 10 memory readings
- °C readings
- Accuracy: $\pm 0.2^{\circ}\text{C}$ from 35.5°C to 42.0°C
- Range of displayed temperature: 34.0°C to 44.0°C
- Operating temperature: 16.0°C to 35.0°C
- Auto shut-off
- Battery life is approximately 4,000 uses
- Supplied with 2*AAA batteries
- 1 year guarantee
- Requires disposable probe covers (see D60.210 / FWH121)

DESCRIPTION	PRODUCT CODE	NPC CODE
Digital Eco Ear Thermometer	D60.200	FWH120
Ear Thermometer Probe Covers	D60.210	FWH121



PROBE COVERS FOR TYMPANIC EAR THERMOMETER ②

- Infection control friendly
- Single-use
- For use with D60.200 & D60.205
- Clip into place for secure fitting
- Pack of 20

DESCRIPTION	PRODUCT CODE	NPC CODE
Ear Thermometer Probe Covers	D60.210	FWH121



SCIAN NON-CONTACT INFRARED THERMOMETER

- Non-contact infrared economy model (LD-FT0100B)
- Non-invasive, infection control friendly
- No probe covers required
- °F & °C readings
- Body & object/ambient temperature modes
- Displayed temperature range:
Body (forehead): 32°C to 43°C
Object/Ambient: 0°C to 70°C
- Accuracy (body mode):
+/- 0.2°C from 35.0°C to 42.0°C
+/- 0.3°C from 32.0°C to 35.0°C & 42.0°C to 43.0°C
- Memory capacity 32 measurements
- Operating temperature range: 10°C to 40°C
- Storage temperature range: -20°C to 55°C
- Supplied with 2x AAA batteries

DESCRIPTION	PRODUCT CODE	NPC CODE
Non-Contact Thermometer	D60.315	FWH384

DIGITAL THERMOMETER ECO/RIGID/FLEXIBLE



- Rigid tip (D60.235 & D60.220)
- Flexible tip (D60.225)
- Easy to read LCD screen
- Measurement response time of 60 seconds (D60.235)
- Measurement response time of 10 seconds (D60.220 & D60.225)
- Compact, accurate & durable
- Oral, axillary & rectal use with bleeper function
- Auto shut-off - 10 minutes
- °C reading (D60.235)
- °F & °C readings (D60.220 & D60.225)
- Supplied with 1.5V button battery
- Size: 134 x 22 x 14mm
- Waterproof (D60.220 & D60.225)

DESCRIPTION	PRODUCT CODE	NPC CODE
Digital Economy Rigid Thermometer	D60.235	FWH259
Digital Rapid Rigid Thermometer	D60.220	FWH224
Digital Rapid Flexible Thermometer	D60.225	FWH123



DIGITAL THERMOMETER COVERS (SHEATHS)



- Universal fitting
- Infection control friendly
- Single-use
- For use with D60.220, D60.225 & D60.235
- Box of 100

DESCRIPTION	PRODUCT CODE	NPC CODE
Digital Thermometer Covers	D60.230	FWH226



EMERALD CLIP-ON ANEROID SPHYGMOMANOMETER

- 2 tube system
- Separate bellows & gauge tubes
- Spring loaded clip-on rear of dial for one-handed use
- 48mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult black nylon calibrated cuff
- Zip-up pouch for convenient storage
- 1 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Emerald Sphygmomanometer	D05.105	FFE183



SAPPHIRE ANEROID SPHYGMOMANOMETER

- 1 tube system
- Combined bellows & gauge tubes
- Easy one-handed use
- 48mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult nylon calibrated cuff
- Zip-up pouch for convenient storage
- 1 year guarantee

COLOUR	PRODUCT CODE	NPC CODE
■ Black	D05.110	FFE184
■ Pink	D05.110P	FFE519
■ Purple	D05.110PL	FFE521



TOPAZ ANEROID SPHYGMOMANOMETER

- 1 tube system
- Combined bellows & gauge tube
- Easy one-handed use
- Trigger button or screw valve pressure release
- Rotating head for left or right handed users
- 54mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult black nylon calibrated cuff
- Zip-up pouch for convenient storage
- 1 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Topaz Sphygmomanometer	D05.115	FFE196

TOPAZ DELUXE ANEROID SPHYGMOMANOMETER



- 1 tube system
- Combined bellows & gauge tube
- Easy one-handed use
- Trigger button or screw valve pressure release
- Robust metal gauge casing
- Rotating head for left or right handed users
- 58mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult black nylon calibrated cuff
- Zip-up pouch for convenient storage
- 2 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Topaz Deluxe Sphygmomanometer	D05.120	FFE522

TOPAZ MULTI-CUFF ANEROID SPHYGMOMANOMETER



- 1 tube system
- Combined bellows & gauge tube
- Easy one-handed use
- Trigger button or screw valve pressure release
- Rotating head for left or right handed users
- 54mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with large adult, adult & child cuffs
- Zip-up pouch for convenient storage
- Supplied in a presentation box
- 1 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Topaz Multi-cuff Sphygmomanometer	D05.130	FFE523

CORAL SHOCK PROOF ANEROID SPHYGMOMANOMETER



- 1 tube system
- Combined bellows & gauge tube
- Easy one-handed use
- Trigger button or screw valve pressure release
- Rotating head for left or right handed users
- Cased within a rubberised shock-proof ring
- Strong nickel-plated metal case
- 59.5mm large dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult black nylon calibrated cuff
- Zip-up pouch for convenient storage
- 2 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Coral Sphygmomanometer	D05.160	FFE819



JADE DESK/WALL ANEROID SPHYGMOMANOMETER

- 2 tube system
- Separate bellows & gauge tubes
- 152mm large dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Desk or wall mounted (supplied with wall bracket)
- Storage for cuff & tubing behind the dial
- 1 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Jade Desk/Wall Sphygmomanometer	D05.150	FFE270



OPAL ANEROID SPHYGMOMANOMETER (FLOOR MODEL)

- 2 tube system
- Separate bellows & gauge tubes
- Stand model
- 147mm large dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Heavy duty plastic feet
- Castors with brakes
- 5 wheels with 360° swivel for ease
- Adjustable height from 95 to 135cm
- Supplied with an adult black nylon calibrated cuff
- Storage basket for cuffs, tubes & in 1 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Opal Sphygmomanometer	D05.170	FFE894



SPHYGMOMANOMETER CUFFS

- Nylon cuffs
- Latex-free
- Push on tube connection (no adaptor required)
- Various sizes from infant to bariatric patients

DESCRIPTION	PRODUCT CODE	NPC CODE
Obese Cuff (40.6 - 66cm), 1 Tube	D06.231	FFE190
Obese Cuff (40.6 - 66cm), 2 Tubes	D06.232	FFE189
Large Adult Cuff (33.3 - 51cm), 1 Tube	D06.201	FFE187
Large Adult Cuff (33.3 - 51cm), 2 Tubes	D06.202	FFE188
Adult Cuff (25.4 - 40.6cm), 1 Tube	D06.221	FFE185
Adult Cuff (25.4 - 40.6cm), 2 Tubes	D06.222	FFE186
Child Cuff (18.4 - 26.7cm), 1 Tube	D06.251	FFE193
Child Cuff (18.4 - 26.7cm), 2 Tubes	D06.252	FFE271
Infant Cuff (13 - 19cm), 1 Tube	D06.261	FFE195
Infant Cuff (13 - 19cm), 2 Tubes	D06.262	FFE194

DIGITAL BLOOD PRESSURE MONITOR



- Fully automatic
- Large LCD display
- One button operation
- IHB Indicator
- British Hypertension Society (BHS) classification
- WHO (blood pressure classification indicator)
- Automatic power off (after 3 minutes of no use)
- Memory capacity of 180 readings
- Powered by 4*AA batteries (supplied)
- 2 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Automatic Digital BP Monitor	D05.205	FFE273
Standard Adult Cuff	D05.206	FFE525
Large Adult Cuff	D05.207	FFE526

DIGITAL DELUX BLOOD PRESSURE MONITOR



- Fully automatic
- Time/date display
- Extra-large LCD display
- British Hypertension Society (BHS) classification
- WHO (blood pressure classification indicator)
- IHB (Irregular heart beat detector)
- Storage compartment behind display
- Ambient temperature display (during sleep mode)
- Low battery indicator
- Automatic power off (after 3 minutes of no use)
- Memory capacity 120 readings
- Powered by 4*AA batteries (supplied)
- 2 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Automatic Digital Delux BP Monitor	D05.230	FFE895
Standard Adult Cuff	D05.206	FFE525
Large Adult Cuff	D05.207	FFE526

DIGITAL WRIST BLOOD PRESSURE MONITOR



- Fully automatic
- Time/date display
- Large LCD display
- British Hypertension Society (BHS) classification
- WHO (blood pressure classification indicator)
- IHB (Irregular heart beat detector)
- Error & low battery indicator
- Automatic power off (after 3 minutes of no use)
- Memory capacity 180 readings
- Powered by 2*AAA batteries (supplied)
- 2 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Automatic Digital Wrist Monitor	D05.240	FFE1207



RUBY SINGLE HEAD STETHOSCOPE

- Quality
- Affordable
- Single head (47mm diameter)
- Ideal for use with a sphygmomanometer cuff
- Aluminium alloy chestpiece
- PVC tube
- 1 year guarantee

COLOUR	PRODUCT CODE	NPC CODE
■ Blue	TD1-13-101	FFE1559
■ Green	TD1-13-104	FFE1560
■ Red	TD1-13-105	FFE1561
■ Black	TD1-13-107	FFE1625
■ Pink	TD1-13-109	N/A



RUBY DUAL HEAD STETHOSCOPE

- Quality
- Affordable
- Dual head chestpiece (47mm diameter)
- Excellent sound transmission
- Lightweight construction
- Aluminium alloy chestpiece
- PVC tube
- 1 year guarantee

COLOUR	PRODUCT CODE	NPC CODE
■ Blue	TD1-13-201	FFE668
■ Green	TD1-13-204	FFE669
■ Red	TD1-13-205	FFE893
■ Black	TD1-13-207	FFE670
■ Purple	TD1-13-208	FFE671
■ Pink	TD1-13-209	FFE672



DUAL HEAD TEACHING STETHOSCOPE

- Dual head chestpiece (47mm diameter)
- Designed specifically for training purposes
- Two binaural headsets coupled with quality tubing
- Aluminium alloy chestpiece
- Black PVC tube
- 1 year guarantee
- Allows the teacher & student to listen to heart, chest & lung sounds simultaneously.

DESCRIPTION	PRODUCT CODE	NPC CODE
Dual Head Teaching Stethoscope	TDI-13-217	FFE673



TWIN TUBE STETHOSCOPE (SPRAGUE RAPPAPORT)

- Dual head chestpiece
- Two tubes connecting the chestpiece to the binaural (45.5mm diameter)
- Provides excellent acoustics in each ear
- Interchangeable chestpieces for adult & paediatrics included
- Zinc alloy chestpiece
- PVC tube
- 1 year guarantee

COLOUR	PRODUCT CODE	NPC CODE
■ Blue	TD1-13-302	FFE674
■ Green	TD1-13-304	FFE675
■ Red	TD1-13-305	FFE676
■ Orange	TD1-13-306	FFE677
■ Black	TD1-13-307	FFE180
■ Purple	TD1-13-308	FFE678
■ Pink	TD1-13-309	FFE679

DIAMOND ADULT STETHOSCOPE

- Dual head chestpiece
- A high quality 48mm diameter chestpiece with a floating diaphragm for excellent acoustics
- Contains internal binaural springs
- Plastic non-chill sleeves for patient comfort
- Supplied with a spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Stainless steel chestpiece
- PVC tube
- 3 year guarantee

COLOUR	PRODUCT CODE	NPC CODE
■ Navy Blue	TD1-13-702	FFE889
■ Burgundy	TD1-13-703	FFE890
■ Hunter Green	TD1-13-704	FFE891
■ Black	TD1-13-707	FFE179





DIAMOND PAEDIATRIC STETHOSCOPE

- Black
- Dual head chestpiece
- Paediatric chestpiece (35mm diameter)
- Floating diaphragm for excellent acoustics
- Contains internal binaural springs
- Plastic non-chill sleeves for patient comfort
- Supplied with spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Stainless steel chestpiece
- PVC tube
- 3 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Diamond Paediatric Stethoscope	TD1-13-717	FFE123



DIAMOND INFANT STETHOSCOPE

- Black
- Dual head chestpiece
- Infant chestpiece (30mm diameter)
- Floating diaphragm for excellent acoustics
- Contains internal binaural springs
- Plastic non-chill sleeves for patient comfort
- Supplied with spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Stainless steel chestpiece
- PVC tube
- 3 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Diamond Infant Stethoscope	TD1-13-727	FFE680



DIAMOND NEONATAL STETHOSCOPE

- Black
- Dual head chestpiece
- Neonatal chestpiece (24mm diameter)
- Floating diaphragm for excellent acoustics
- Contains internal binaural springs
- Plastic non-chill sleeves for patient comfort
- Supplied with a spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Stainless steel chestpiece
- 3 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Diamond Neonatal Stethoscope	TD1-13-737	FFE122

CARDIOLOGY ADULT STETHOSCOPE



- Dual head chestpiece
- 48mm diameter chestpiece
- Deep cone-shaped bell for excellent acoustic sensitivity
- High quality tube for resistance to ambient interference
- A-shaped tube design features a binaural with reinforced stainless steel inner-spring for durability
- Features an anti-chill rim for patient comfort
- Supplied with a spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Stainless steel chestpiece
- PVC tube
- 3 year guarantee

COLOUR	PRODUCT CODE	NPC CODE
■ Navy Blue	TD1-13-902	FFE1210
■ Black	TD1-13-907	FFE820

DELUXE CARDIOLOGY ADULT STETHOSCOPE



- Ergonomic & unique design
- 50mm diameter chestpiece
- Quality tube for excellent acoustic performance
- Comfortable & flexible binaural
- Durable, soft silicone ear tips
- Supplied with a spare diaphragm & spare ear tips
- Supplied in a foam lined box
- Zinc alloy chestpiece
- PVC tube
- 3 year guarantee

COLOUR	PRODUCT CODE	NPC CODE
■ Black	TD1-13-807	FFE681

PINARD STETHOSCOPE



- For listening to fetal heart sounds
- Available in plastic, aluminium or wood
- Approximately 14.5cm

COLOUR & SIZE	PRODUCT CODE	NPC CODE
■ Plastic Black	03.350.10	FFE682
■ Plastic Red	03.350.15	FFE683
■ Plastic Pink	03.350.20	FFE684
■ Aluminium	03.360.20	FFE685
■ Wood	03.370.20	FFE686



EMERALD SPHYGMOMANOMETER & SPRAGGUE STETHOSCOPE SET

Sphygmomanometer

- Separate bellows & gauge tubes
- Spring loaded clip-on rear of dial for one-handed use
- 48mm dial for easy reading of pressure values
- Push on tube connection (no adaptor required)
- Supplied with an adult nylon calibrated cuff

Stethoscope

- Dual head chestpiece
- 2 tubes connecting the chestpiece to the binaural
- Interchangeable chestpieces for adult & paediatrics incl.

DESCRIPTION	PRODUCT CODE	NPC CODE
Black Emerald Aneroid Sphygmomanometer & Black Twin Tube Sprague Rappaport Stethoscope Set	D05.305	FFE823
Pink Emerald Aneroid Sphygmomanometer & Pink Twin Tube Sprague Rappaport Stethoscope Set	D05.305P	FFE824
Purple Emerald Aneroid Sphygmomanometer & Purple Twin Tube Sprague Rappaport Stethoscope Set	D05.305PL	FFE825

MIDWIFERY SET PINK (COMPLETE)

Includes:

- Stethoscope
- Sphygmomanometer
- Pinard
- Tourniquet
- Fob watch
- Pen torch
- Face shield
- Dressing scissors
- Hand sanitiser
- Utility scissors
- Obstetric wheel



DESCRIPTION	PRODUCT CODE	NPC CODE
Midwifery Set Pink (Complete)	D05.905	FAK112
Pink Plastic Case	D05.950	N/A

NURSES SET PINK (COMPLETE)

Includes:

- Stethoscope
- Sphygmomanometer
- Tourniquet
- Fob watch
- Pen torch
- Face shield
- Dressing scissors
- Hand sanitiser
- Utility scissors



DESCRIPTION	PRODUCT CODE	NPC CODE
Nurses Set Pink (Complete)	D05.910	FAK111
Pink Plastic Case	D05.950	N/A

PATELLA HAMMER



02.290.40

- Reflex Hammers used to test deep tendon reflexes
- Perfect for neurological physical examinations
- Available in a variety of patterns

DESCRIPTION	PRODUCT CODE	NPC CODE
Percussor with Telescopic Handle	02.270.50	FFE1793
Buck Percussor (19cms) Zinc Alloy Handle	02.290.40	FFE1794
Taylor Percussor S/Steel (19cms) Zinc Alloy Handle	02.340.20	FFE1795
Babinski Percussor (23cms) Zinc Alloy Handle	02.400.20	FFE1796
NHP (Queens Square) Percussor (32cm) Nylon Handle	02.450.20	FBF2351



02.270.50



02.340.20



02.400.20



02.450.20

GARDINER BROWN TUNING FORK



- Classic Gardiner Brown design tuning fork
- High quality surgical steel tuning fork medical use
- Base plate for Rinne and Weber testing in medicine
- 6 base frequencies available

DESCRIPTION	PRODUCT CODE	NPC CODE
Tuning Fork 128hz C/S	02.560.10	FFC800
Tuning Fork 256hz C/S	02.560.20	FFC803
Tuning Fork 512hz C/S	02.560.30	FFC804
Tuning Fork 1024hz C/S	02.560.40	FFC798

RING CUTTER AND BLADES



- The high quality Timesco ring cutters
- Supplied with spare blades
- Essential medical tool for the cutting/removal of rings
- Made from high grade stainless steel

DESCRIPTION	PRODUCT CODE	NPC CODE
Ring Cutter & Blade	02.670.20	FFQ5256
Ring Cutter Blade Pack of 10	02.670.50	N/A



Bayonet fitting



Pin fitting

ORION XENON STANDARD OPHTHALMOSCOPE

- Available in pin (screw) or bayonet fitting
- Lightweight hermetically sealed durable plastic head
- Precision optics
- Latex-free rubber eyebrow rest
- Bright white light from long life Xenon / standard bulbs
- Clear illumination in ophthalmic diagnostic procedures
- Supplied in a zip-up case
- 3 year guarantee (excludes consumables)

DESCRIPTION	PRODUCT CODE	NPC CODE
Ophthalmoscope - Xenon (pin)	2501.250.50	FFC146
Ophthalmoscope - Xenon (bayonet)	2501.400.50	FER14211
Spare bulb, No. 80, 2.5v	Z-2700.080.21	FFC144



Bayonet fitting



Pin fitting

ORION XENON OTOSCOPE

- Available in pin (screw) or bayonet fitting
- Durable metal handle & head
- Non-reflective matt finish
- Bright white light from long life Xenon bulbs
- Universal fittings for aural & nasal
- Supplied in a zip-up case
- 3 year guarantee (excludes consumables)

DESCRIPTION	PRODUCT CODE	NPC CODE
Otoscope - Xenon (pin)	2500.100.50	IVN012
Otoscope - Xenon (bayonet)	2500.400.50	FER14210
Spare bulb, No. 50, 2.5v	Z-2700.050.21	FFC143
Adaptor for disposable specula	2500.810.50	N/A



SINGLE-USE EAR SPECULA



- Single-use
- Available in 2.5mm & 4mm
- Available in packs of 100 or bulk carton of 10,000

DESCRIPTION	PRODUCT CODE	NPC CODE
Single-use specula - 2.5mm (x 100)	7000.850.55	IVN003
Single-use specula - 2.5mm (x 10,000)	7000.850.25	FER14205
Single-use specula - 4mm (x 100)	7000.850.60	FBF2348
Single-use specula - 4mm (x 10,000)	7000.850.40	FER14206

ORION DIAGNOSTIC SET



- Available in pin (screw) or bayonet fitting
- 3 year guarantee (excludes consumables)
- Supplied in hard plastic case

Contents:

- Interchangeable handle
- Otoscope head
- Ophthalmoscope head
- Spare otoscope bulb
- Spare ophthalmoscope bulb
- Reusable specula sizes 1, 2 & 3

DESCRIPTION	PRODUCT CODE	NPC CODE
Set - Xenon (pin)	2502.275.75	IVN013
Set - Xenon / standard (bayonet)	2502.400.50	FER14212

ORION XENON DELUXE SET



- Available in pin (screw) fitting only
- 3 year guarantee (excludes consumables)
- Supplied in hard plastic case

Contents:

- Interchangeable handle
- Otoscope head
- Ophthalmoscope head
- Disposable speculum
- Disposable speculum adaptor
- Reusable specula 2.5mm, 3.5mm, 4.5mm
- Nasal speculum screw
- Spare otoscope bulb
- Spare ophthalmoscope bulb

DESCRIPTION	PRODUCT CODE	NPC CODE
Deluxe set - Xenon (pin)	2502.275.90	IVN014

ORION XENON COMBO SET



- Available in pin (screw) fitting only
- 3 year guarantee (excludes consumables & case)
- Supplied in zip-up case

Contents:

- Interchangeable handle
- Otoscope head
- Ophthalmoscope head
- Reusable specula 2.5mm, 3.5mm, 4.5mm
- Nasal speculum screw
- Laryngeal stem
- Mirror No.3
- Mirror No.4
- Depressor holder
- Spare otoscope bulb

DESCRIPTION	PRODUCT CODE	NPC CODE
Combo set - Xenon (pin)	2502.300.50	FFC802



OPTIMA MINI POCKET OTOSCOPE

- Available in Xenon or LED versions
- Durable & lightweight
- Supplied in a soft pouch
- A fibre optic light source
- Bright white light from long life bulbs
- 3 year guarantee (excludes consumables)

DESCRIPTION	PRODUCT CODE	NPC CODE
Otoscope - Xenon, with spare bulb	2505.100.25	IVN022
Otoscope - Xenon	2505.105.25	FER14202
Spare bulb - Xenon	2505.325.11	FFC142
Otoscope - LED	2505.500.50	FER14200
Disposable specula - 2.5mm (x 100)	7000.850.55	IVN003
Disposable specula - 4mm (x 100)	7000.850.60	FBF2348



OPTIMA MINI POCKET OPHTHALMOSCOPE

- Durable & lightweight
- Supplied in a soft pouch
- Long life bright white Xenon bulb
- Clear illumination in ophthalmic diagnostic procedures
- Macular beam dioptries lenses from 0 to +20 & 0 to -20
- 3 year guarantee (excludes consumables)

DESCRIPTION	PRODUCT CODE	NPC CODE
Ophthalmoscope, with spare bulb	2505.250.25	FFE1631
Ophthalmoscope	2505.255.25	FER14203
Spare bulb	2505.325.21	IVN023



OPTIMA MINI POCKET DIAGNOSTIC SET

- Available in Xenon or LED versions (otoscope only)
- 3 year guarantee (excludes consumables)
- Supplied in hard plastic case

Contents:

- Ophthalmoscope complete with handle (Xenon)
- Otoscope complete with handle (LED or Xenon)
- Includes disposable specula (7 x 2.5mm / 7 x 4.0mm)

DESCRIPTION	PRODUCT CODE	NPC CODE
Pocket set - Xenon, with spare bulbs	2505.275.75	IVN021
Pocket set - Xenon	2505.275.50	FER14204
Pocket set - LED/Xenon	2505.575.75	FER14201



OPTIMA NEO OTOSCOPE

- Superb bright Xenon white light illumination
- Can be used with disposable & reusable specula which are available in various diameters
- A 3x magnification swivel lens is incorporated in the otoscope head which when rotated allows it to be used for aural examinations
- The otoscope head can also be used as an oral illuminator
- An insufflation port is also provided in the otoscope head for connection to insufflation bellows for pneumatic testing of the tympanic membrane
- Supplied in a hard case with a spare bulb & a set of reusable specula & disposable specula
- 2 year guarantee (excludes consumables)

DESCRIPTION	PRODUCT CODE	NPC CODE
Otoscope (2.5v)	3600.100.10	FER14207
Spare bulb	3700.100.10	N/A



Apertures



OPTIMA NEO OPHTHALMOSCOPE

- Superb bright Xenon white light illumination
- Various aperture sizes
- Apertures:
 - Micro spot: for very small undilated pupils
 - Small spot: for undilated pupils, ideal aperture to start visuals
 - Standard spot: for dilated pupils, general application
 - Fixation aperture: for observation of eccentric fixation
 - Slit aperture: for observation of lesions
 - Red free aperture: for observation of retinal abnormalities
- Supplied in a hard case with a spare bulb
- 2 year guarantee
- There are 28 colour coded lenses, green for + & red for -

0	1	2	3	4	5	6	7	8	9	10	12	15	20	40
	1	2	3	4	5	6	7	8	9	10	12	20	25	

DESCRIPTION	PRODUCT CODE	NPC CODE
Ophthalmoscope (2.5v)	3600.110.10	FER14208
Spare bulb	3700.110.10	N/A



OPTIMA NEO DIAGNOSTIC SET

- Contents:
- Optima Neo ophthalmoscope complete with handle
 - Optima Neo otoscope complete with handle
 - Reusable specula (2.5 / 3.5 / 4.5 / 5.5 / 9.0mm)
 - Single-use specula (7 x 2.5mm & 7 x 4.0mm)
 - Spare bulbs for ophthalmoscope & otoscope
 - 2 year guarantee (excludes consumables)

DESCRIPTION	PRODUCT CODE	NPC CODE
Diagnostic set (2.5v)	3600.120.10	FER14209



DISPOSABLE PEN TORCH



- Single-use
- Pocket clip
- Plastic sealed casing
- Non-replaceable battery
- Features a pupil gauge
- Integrated pocket clip, when pressed powers the torch
- Supplied in packs of six

DESCRIPTION	PRODUCT CODE	NPC CODE
Disposable pen torch	D85.100	FFE1811



PLASTIC REUSABLE PEN TORCH

- Reusable (batteries can be replaced)
- Plastic casing
- Integrated pocket clip, when pressed powers the torch
- Supplied with 2*AAA batteries
- Provided in a blister presentation pack

COLOUR	PRODUCT CODE	NPC CODE
■ Blue	D85.005	FFE1222
■ Green	D85.010	FFE1223
■ Red	D85.015	FFE1224
□ White	D85.020	FBF2352
■ Pink	D85.025	FFE1225
■ Purple	D85.030	FFE1226



DELUXE STAINLESS STEEL REUSABLE PEN TORCH

- Reusable (batteries can be replaced)
- Stainless steel casing
- Features a pupil gauge
- Pocket clip
- Supplied with 2*AAA batteries
- Provided in a blister presentation pack

DESCRIPTION	PRODUCT CODE	NPC CODE
Deluxe pen torch	D85.300	FFE1227



EASY TO USE & TRACK

Intuitive controls for easy to understand operation. Automated QPI score measuring patient pupillary reflexes. Track trends and changes to detect early warning signs.



IDENTIFY NEUROLOGICAL ANOMALIES

Track precise pupil measurements (miosis/mydriasis). Aid in early detection of neurological anomalies such as anisocoria.

PORTABLE & NON-INVASIVE

Portable hand-held design for easy use throughout the facility. Wipeable, re-usable silicon eye cup for maximum patient comfort.



NEUROLIGHT PUPILLOMETER

- Portable Design
- Easy to use interface
- Opaque reusable eye cup
- Accurately track pupil diameter and response time
- Output measurements into easy to read QPI score
- Identify neurological anomalies via pupil
- Log and track patient specific data
- Transferable data for further analysis

Neurolight Pupilometer / D85.500
supplied with:

- Neurolight
- Wireless charging station
- 10 eyecups
- Power supply

DESCRIPTION	PRODUCT CODE	NPC CODE
Neurolight Pupilometer	D85.500	FBF85057
Neurolight Without Accessories	D85.505	N/A
Neurolight Eye Cups Box of 10	D85.510	N/A



UTILITY SCISSORS

- Tuff cuts
- Stainless steel blades
- Tough/durable plastic handle
- Available in a variety of colours
- Lower blade has a serrated edge & is hooked & rounded at the lip for ease of use
- Available in two sizes, 19cm/7.5" & 15cm/6"
- Pack of 10

COLOUR & DESCRIPTION	PRODUCT CODE	NPC CODE
Utility Scissors 7.5" Green	12.711.06	FFQ4563
Utility Scissors 7.5" Yellow	12.711.11	FFQ4608
Utility Scissors 7.5" Orange	12.711.16	N/A
Utility Scissors 7.5" Blue	12.711.21	FFQ4564
Utility Scissors 7.5" Black	12.711.31	FFQ4565
Utility Scissors 7.5" Red	12.711.36	FFQ4566
Utility Scissors 7.5" Pink	12.711.41	FFQ4605
Utility Scissors 7.5" Purple	12.711.46	FFQ4606

COLOUR & DESCRIPTION	PRODUCT CODE	NPC CODE
Utility Scissors 6" Green	12.751.06	FDB85059
Utility Scissors 6" Blue	12.751.21	FDB85060
Utility Scissors 6" Black	12.751.31	FDB85061
Utility Scissors 6" Red	12.751.36	FDB85062
Utility Scissors 6" Pink	12.751.41	FDB85063
Utility Scissors 6" Purple	12.751.46	FDB85064



Pin clasp

Battery save clip

SILICONE FOB WATCH

- New improved design
- Japanese Quartz movement
- Pin clasp closure
- Made from silicone based material
- Ideal for infection control
- Contains 1 x watch battery (SR626SW)
- 2 year battery guarantee

COLOUR	PRODUCT CODE	NPC CODE
■ Pink	TDM-FW-500	FFE1214
■ Blue	TDM-FW-505	FFE1209
■ Green	TDM-FW-510	FFE1215
■ Red	TDM-FW-515	FFE1216
■ White	TDM-FW-520	FFE1217
■ Black	TDM-FW-525	FFE1218
■ Purple	TDM-FW-530	FFE1219
■ Pack of 14 (2 of each colour)	TDM-FW-545	FFE1220



CARABINER WATCH

- Carabiner clip-on style
- Robust metal casing
- Clip to prevent battery engagement prior to use
- Japanese Quartz movement
- Contains 1 x watch battery (SR626SW)
- 2 year battery guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Carabiner fob watch	TDM-FW-550	FFE1221



WOODEN TONGUE DEPRESSOR



- Strong & Flexible
- Smooth rounded tips
- Dimensions: 150mm x17mm
- Non-Sterile
- FSC Certified

DESCRIPTION	PRODUCT CODE	NPC CODE
Wooden Tongue Depressor (Pack of 100)	D15.100	FBF85056



STANDARD TOURNIBUTTON TOURNIQUET



- Single-use to reduce the risk of cross infection
- Supplied in a convenient dispensing box
- Easy to reposition & tighten with a button/stud fastener
- Adult & Paediatric use
- Material: TPE band
- Latex-free
- Band size: 381mm x 25mm
- Box of 100

DESCRIPTION	PRODUCT CODE	NPC CODE
Standard TourniButton Single-Use Tourniquets (Button Type)	D80.200	FWJ103

NEW
PRODUCT

PAEDIATRIC TOURNIBUTTON TOURNIQUET



- Single-use to reduce the risk of cross infection
- Supplied in a convenient dispensing box
- Easy to reposition & tighten with a button/stud fastener
- Paediatric use
- Material: TPE band
- Latex-free
- Band size: 255mm x 18mm
- Box of 200

DESCRIPTION	PRODUCT CODE	NPC CODE
Paediatric TourniButton Single-Use Tourniquets (Button Type)	D80.205	FVH85097

NEW
PRODUCT

BARIATRIC TOURNIBUTTON TOURNIQUET



- Single-use to reduce the risk of cross infection
- Supplied in a convenient dispensing box
- Easy to reposition & tighten with a button/stud fastener
- Bariatric use
- Material: TPE band
- Latex-free
- Band size: 812mm x 35mm
- Box of 30

DESCRIPTION	PRODUCT CODE	NPC CODE
Bariatric TourniButton Single-Use Tourniquets (Button Type)	D80.210	FVH85096



TOURNIBAND TOURNIQUET



- Single-use to reduce the risk of cross infection
- Supplied on a perforated roll
- Stretch band
- Supplied in dispenser pack
- Latex-free
- Material: TPE
- Band size: 45cm x 2.5cm (approx.)
- Available in boxes of 25 or 100

DESCRIPTION	PRODUCT CODE	NPC CODE
TourniBand Single Use Tourniquets Box 25	D80.100	FWJ195
TourniBand Single Use Tourniquets Box 100 (Pre-cut, Not on Roll)	D80.101	FWJ196



TOURNICLIP QUICK RELEASE TOURNIQUET

- Reusable
- Non pinch for patient comfort
- Quick clip application & button release for ease of use
- Available in a variety of colours
- Latex-free

COLOUR	PRODUCT CODE	NPC CODE
■ Blue	D80.005	FWJ203
■ Green	D80.010	FWJ205
■ Red	D80.015	FWJ207
■ Pink	D80.020	FWJ202
■ Purple	D80.025	FWJ206
■ Yellow	D80.030	FWJ204



FINGER TIP PULSE OXIMETER CN130

- Suitable for adult & paediatric use (2 years+)
- 30 hours continuous use
- One colour LED display
- Displays SpO2, pulse rate & pulse bar
- Supplied with batteries (2*AAA), protective rubber boot & lanyard
- Pulse rate measurement range:
30-99bpm: ± 2 bpm accuracy / 100-250bpm: $\pm 2\%$ accuracy
- SpO2 measurement range: 70% - 100%: $\pm 2\%$ accuracy
- 1 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
CN130 Pulse Oximeter	TDM-PO-CN130	FBC1668
Pulse Oximeter Carry Case	TDM-PO-C3CC	N/A



FINGER TIP PULSE OXIMETER CN330

- Suitable for adult & paediatric use (2 years+)
- 30 hours continuous use
- Dual colour OLED screen, displays SpO2, pulse rate & Plethysmogram waveform
- 4 directional displays
- Supplied with 2*AA batteries
- Pulse rate measurement range:
30-99bpm: ± 2 bpm / 100-250bpm: $\pm 2\%$
- SpO2 measurement range:
70% - 100%: $\pm 2\%$ accuracy
- 1 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
CN330 Pulse Oximeter	TDM-PO-CN330	FBC85144
Pulse Oximeter Carry Case	TDM-PO-C3CC	N/A



FINGER TIP PULSE OXIMETER C19

- Compact & durable
- Suitable for adult & paediatric use (2 years+)
- 30 hours continuous use
- One colour LED display
- Displays SpO2, pulse rate & pulse bar
- Supplied with batteries (2*AAA), carry case, protective rubber boot & lanyard
- Pulse rate measurement range:
30-99bpm: ± 2 bpm accuracy / 100-235bpm: $\pm 2\%$ accuracy
- SpO2 measurement range: 70% - 99%: $\pm 3\%$ accuracy
- 2 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
C19 Pulse Oximeter	TDM-PO-C19	FBC323
Pulse Oximeter Carry Case	TDM-PO-C3CC	N/A



FINGER TIP PULSE OXIMETER C29

- Suitable for adult & paediatric use (2 years+)
- 30 hours continuous use
- Dual colour OLED screen displays SpO2, pulse rate, pulse bar & plethysmogram waveform
- Adjustable brightness & 6 display modes
- Auto power down & low battery indicator
- Supplied with batteries (2*AAA), carry case, protective rubber boot & lanyard
- Pulse rate measurement range:
30-99bpm: ± 2 bpm accuracy / 100-235bpm: ± 2 % accuracy
- SpO2 measurement range: 70% - 100%: ± 2 % accuracy
- 2 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
C29 Pulse Oximeter	TDM-PO-C29	FBC385
Pulse Oximeter Carry Case	TDM-PO-C3CC	N/A



PAEDIATRIC FINGER TIP PULSE OXIMETER C52

- Small finger placement area for an accurate reading
- For paediatrics & adults (2+ years), compact for accuracy
- 30-hour use, auto power down, low battery indicator
- Dual-colour OLED: SpO2, pulse rate, pulse bar, waveform
- Adjustable brightness, 6 display modes
- Includes 2x AAA batteries, carry case, lanyard
- Pulse Rate: 30-99bpm (± 2 bpm), 100-235bpm (± 2 %) | SpO2: 70%-100% (± 2 %)
- 3-year warranty

DESCRIPTION	PRODUCT CODE	NPC CODE
C52 Pulse Oximeter	TDM-PO-C52	FBC269
Pulse Oximeter Carry Case	TDM-PO-C3CC	N/A



FINGER TIP PULSE OXIMETER CI218

- For adults & children (2+ years)
- Durability: Drop-proof (1m), waterproof (level 2)
- Display: OLED shows SpO2, pulse rate, pulse bar & perfusion index
- Power: Auto power down, low battery indicator, battery status, 35-hour continuous use, low power consumption
- Features: Weak signal prompt, data storage
- Includes: 2x AAA batteries, deluxe case & lanyard
- Measurement Ranges & Accuracy:
 - Pulse Rate: 30-99bpm (± 2 bpm), 100-250bpm (± 2 %)
 - SpO2: 70%-100% (± 2 %)
 - Perfusion Index: 0.1-20% (0.1% resolution)
 - Low Perfusion Accuracy: SpO2 70%-100% (± 3 %), PR 30-250bpm (± 2 %)
- Warranty: 5 years

DESCRIPTION	PRODUCT CODE	NPC CODE
CI218 Pulse Oximeter	TDM-PO-CI218	FBC1667



WRIST PULSE OXIMETER

- Real time and sleep measurement modes
- Spot Check reading / 20 hour continuous reading
- OLED touch screen
- Bluetooth / transfer data to smartphone app
- Stores 72 hours of sleep data
- Lithium-ion rechargeable battery
- Pulse rate measurement range:
30-99bpm: ± 2 bpm accuracy / 100-250bpm: $\pm 2\%$ accuracy
- SpO2 measurement range:
70% - 100%: $\pm 2\%$ accuracy
- 1 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Wrist Pulse Oximeter	TDM-PO-300W628	FBC85143



300M HANDHELD PULSE OXIMETER

- Suitable for adult, paediatric & neonatal use
- 16 hours of continuous use
- High resolution 2.8" TFT colour display
- Displays SpO2, pulse rate, pulse bar & plethysmogram
- Adjustable audio & visual alarms
- Patient ID ranges from 1 to 127
- Auto power off & low battery indicator
- Supplied with batteries (3*AA)
- Pulse rate measurement range:
30-235bpm: ± 2 bpm or 2% accuracy (the largest)
- SpO2 measurement range: 70% - 100%: $\pm 2\%$ accuracy
- 3 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
300M Handheld Pulse Oximeter	TDM-PO-300M	FBC270





Adult



Paediatric



Neonatal

FINGER CLIP SP02 SENSORS (FOR 300M)

- Reusable
 - Available for adult use 30kg>
- Available for paediatric use 10-50kg
 - ‘Y’ type available for neonatal use <3kg

DESCRIPTION	PRODUCT CODE	NPC CODE
Adult Finger SP02 Sensor Probe	TDM-PO-305	FBC307
Paediatric Finger SP02 Sensor Probe	TDM-PO-306	FBC308
Neonatal Finger SP02 Sensor Probe	TDM-PO-307	FBC309



Adult



Infant



Paediatric

SINGLE-USE MICRO FOAM FINGER SENSORS (FOR 300M)



- Single-use
 - Finger probes for the 300M hand held unit
- Available in adult, paediatric & neonate sizes

DESCRIPTION	PRODUCT CODE	NPC CODE
Adult Single-Use Micro Foam Finger Sensor	TDM-PO-325	FBC1664
Paediatric Single-Use Micro Foam Finger Sensor	TDM-PO-326	FBC1665
Infant Single-Use Micro Foam Finger Sensor	TDM-PO-327	FBC1666



★ NEW
PRODUCT

HI-BEBE BT-200L FOETAL DOPPLER

- Compact & lightweight
- Foetal heart rate processing
- High resolution LCD display
- Built in loudspeaker
- 2 Mhz frequency
- Suitable for water birth delivery - IPX7 waterproof probe
- Auto power off
- Low battery indicator
- Battery life: 5 hours continuous use
- Pulse indicator
- Measuring range: 50bpm to 240bpm
- Weight: 190g without batteries
- FHR & sound management
- Supplied with carry case
- 1 Year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Hi-Bebe 200L Foetal Doppler	TDM-BT-200L	FGR1023
Hi-Bebe Mobile Doppler Stand	TDM-BT-230	N/A
Hi-Bebe 200L Spare Doppler Probe	TDM-BT-200LP	N/A



★ NEW
PRODUCT

HI-BEBE BT-220C FOETAL DOPPLER

- Compact & lightweight
- Foetal heart rate processing
- 2.4" Colour LCD display with FHR trend display
- Built in loudspeaker
- 2 MHz frequency
- Suitable for water birth delivery - IPX7 waterproof probe
- Auto power off
- Low battery indicator
- Battery life: 4 hours continuous use
- Pulse indicator
- Measuring range: 50bpm to 240bpm
- Weight: 200g without batteries
- Acoustic stimulator function
- BMI & body fat analysis
- FHR & sound management
- Supplied with carry case
- 1 year guarantee

DESCRIPTION	PRODUCT CODE	NPC CODE
Hi-Bebe 220C Foetal Doppler	TDM-BT-220C	FGR1024
Hi-Bebe Mobile Doppler Stand	TDM-BT-230	N/A
Hi-Bebe 220C Spare Doppler Probe	TDM-BT-220CP	N/A



SECA 384 BABY & TODDLER DIGITAL WEIGHING SCALES

- Capacity: 20kg
- Graduation: 10 g < 10 kg > 20 g
- Dimensions: 552 x 162 x 362mm
- Weight: 2.8kg
- Power supply: 6*AA batteries
- Functions: TARE / weighing range switch-over / automatic switch-off / BMIF / hold

DESCRIPTION	PRODUCT CODE	NPC CODE
Seca 384 Baby & Toddler Scales	TDM-S-384	N/A



Illustrated without tray

SECA 385 BABY, TODDLER & CHILD DIGITAL WEIGHING SCALES

- Capacity: 50kg
- Graduation: 20 g < 20 kg > 50 g
- Dimensions: 552 x 162 x 362mm
- Weight: 2.8kg
- Power supply: 6*AA batteries
- Functions: TARE / weighing range switch-over / automatic switch-off / BMIF / hold

DESCRIPTION	PRODUCT CODE	NPC CODE
Seca 385 Baby, Toddler & Child Scales	TDM-S-385	N/A



SECA CARRY CASE FOR 384 & 385 SCALES

- For transporting Seca 384 & 385 scales.
- Provides optimum protection
- Large opening
- Adjustable shoulder strap
- Water-repellent nylon surface

DESCRIPTION	PRODUCT CODE	NPC CODE
Seca Carry Case for 384/385	TDM-S-413	N/A



SECA BACK PACK FOR 384 & 385 SCALES

- For transporting Seca 384 & 385 scales.
- Provides optimum protection
- Large opening
- Adjustable shoulder straps
- Water-repellent nylon surface

DESCRIPTION	PRODUCT CODE	NPC CODE
Seca Backpack for 384/385	TDM-S-413B	N/A



COMPRESSOR NEBULISER NB-210C

- Compact & streamlined design
- Built-in handle for easy portability
- Ideal for home, clinical or hospital use
- Features a convenient storage compartment for storing the chamber, masks, tubing & other spares
- Mains operated
- Medicine chamber maximum capacity: 6ml
- Particle size: MMAD 5 microns
- Nebulisation rate: $\geq 0.2\text{ml} / \text{min}$
- Noise level: $\leq 65\text{dB}$
- Dimensions: 330mm x 200mm x 138mm

Includes:

Adult mask, Child mask, Mouthpiece, Filter cotton x 5, Connecting tubing, Nebulising chamber

DESCRIPTION	PRODUCT CODE	NPC CODE
Scian Compressor Nebuliser NB-210C	D10.000	FDD4178
Nebuliser Filters (Pack of 5)	D10.001	FDK85415



COMPRESSOR NEBULISER NB-216C

- Pressure Adjustable
- LED Display
- Time Countdown
- Low Residue
- Size: 166*136*104mm
- Power Input: DC 12.0V, 1.35A
- Pressure: 30-50psi
- Free Flow Range: $>4\text{L/min}$
- Nebulising Pressure: $>6\text{psi}$
- Noise Level (with Nebuliser Cup): $< 55\text{dB}$
- Operation Mode: Continuous Work

DESCRIPTION	PRODUCT CODE	NPC CODE
Scian Compressor Nebuliser NB-216C	D10.002	N/A
Nebuliser Filters (Pack of 5)	D10.001	FDK85415

★ NEW PRODUCT

NEBULISER KIT WITH MOUTHPIECE

2



- Supplied with; O2 tubing & nebuliser medication chamber (6ml & mouthpiece)
- Light weight medical grade PVC
- Latex & Phthalates free
- Box of 50

DESCRIPTION	PRODUCT CODE	NPC CODE
Nebuliser Kit with Mouthpiece	D10.017	FAG85148

★ NEW
PRODUCT

NEBULISER KIT

| 2 |

- Supplied with; mask, O2 tubing & nebuliser medication chamber (6ml)
- With nose clip & elastic head strap
- Light weight medical grade PVC
- Latex & Phthalates free
- Box of 50

DESCRIPTION	PRODUCT CODE	NPC CODE
Nebuliser Kit Adult	D10.015	FAG7011
Nebuliser Kit Child	D10.016	FAG7018

★ NEW
PRODUCT

NON-REBREATH MASK

| 2 |

- Supplied with; mask, O2 tubing & reservoir bag
- With nose clip & elastic head strap
- Light weight medical grade PVC
- Latex & Phthalates free
- High concentration 99% to 100% delivery
- Box of 50

DESCRIPTION	PRODUCT CODE	NPC CODE
Adult Non-Rebreath Mask	D10.305	FAG7028
Paediatric Non-Rebreath Mask	D10.306	FAG7025

★ NEW
PRODUCT

MEDIUM CONCENTRATION FACE MASKS

| 2 |

- Supplied with; mask & O2 tubing
- With nose clip & elastic head strap
- Light weight medical grade PVC
- Latex & Phthalates free
- Medium concentration 35% to 50% delivery
- Box of 50


DESCRIPTION	PRODUCT CODE	NPC CODE
Adult Medium Concentration Mask	D10.315	FAG7031
Paediatric Medium Concentration Mask	D10.316	FDD85369



EVAC-S SINGLE-USE SUCTION DEVICE

| 2 |

- Manual handheld suction device
- Completely single-use
- Easy to assemble
- Product is ready to use after connecting the catheter
- After use the device can be placed into the polybag, resealed & the product label can be peeled to reveal a biohazard warning label
- Vacuum (max): 500mmHg
- Peak free airflow: >20ml
 - Disposable container volume: 300ml
 - External diameter of cap (suction port): 17mm
 - External diameter of catheters:
 - 9mm Adult, 10fg Paediatric
- Operating environmental temperature: -20°C to +50°C
- Storage environmental temperature: -40°C to +60°C
- Net weight: 270g



- The box contents are:
- Suction pump handle
 - 300ml canister
 - Adult catheter
 - Paediatric catheter
 - Self-seal disposal polybag
 - Bio-hazard label ready to be displayed for disposal
 - Multilingual product instructions
 - Multilingual product labeling
 - Supplied in boxes of 5

DESCRIPTION	PRODUCT CODE	NPC CODE
EVAC-S suction device	TDM-ESD-100	FPA626



RESCUER MVP SUCTION DEVICE

| 2 |

- Reusable handle
- Single-use canister(s) / catheter(s)
- Lightweight, portable emergency device
- Easy one-handed use
- 100° swivel allowing the operator greater accessibility in tight situations & avoiding premature engagement of the overflow mechanism
- Provides efficient operation & simple maintenance
- 50% to 100% suction setting to allow the device to be used on neonates to large adults

DESCRIPTION	PRODUCT CODE	NPC CODE
Rescuer MVP Suction Device	TDM-MVP-100	FPA635
Spare Collection Jar / Catheters	TDM-MVP-200	FPA636

FACE SHIELD IN KEYRING POUCH



- Single-use
- Available in various coloured keyring pouches
- Convenient & effective protection from cross-infection during mouth-to-mouth resuscitation
- One-way valve filter

COLOUR	PRODUCT CODE	NPC CODE
Blue	TDM-PM-051	FDG1622
Black	TDM-PM-052	FDG1623
Purple	TDM-PM-053	FDG1624
Red	TDM-PM-054	FDG1625
Yellow	TDM-PM-056	FDG1626
Orange	TDM-PM-057	FDG1627
Green	TDM-PM-058	FDG1628
Pink	TDM-PM-059	FDG1629
Pack of 40 (5 of each colour)	TDM-PM-045	N/A



TIMESCO DISPOSABLE MANIKIN FACE SHIELDS



- Disposable manikin face shields for use in mouth-to-mouth resuscitation/CPR training situations
- Helps ensure student safety
- Supplied in rolls of 50

DESCRIPTION	PRODUCT CODE	NPC CODE
Timesco Disposable Manikin Face Shields Roll of 50	TDM-PM-175	FAG85223



FACE SHIELD IN SACHET



- Face shield for personal protection during mouth-to-mouth resuscitation
- Supplied in easy tear sachet
- Paper filter

DESCRIPTION	PRODUCT CODE	NPC CODE
Face shield in sachet	TDM-PM-165	FDG872



CPR RESCU-MASK

2

- Single-use
- Replacement filters available for training purposes

Supplied with:

- One-way bacterial & viral filter
- Built-in O2 port
- Elasticated head strap
- Hard plastic carry case

DESCRIPTION	PRODUCT CODE	NPC CODE
CPR Rescu-Mask	TDM-PM-250	FDD2552
One-way valve / filter (x5)	TDM-PM-110	FDA85091



SINGLE-USE RESUSCITATORS (BVM)

2

- Textured grip (non-slip)
- 360° swivel connector for greater versatility
- Glued components (with exception of face mask)
- Adjustable air cushion face mask for maximum patient comfort
- Transparent bag for complete internal visibility
- Excellent recoil recovery following bag compression
- Available in a variety of sizes
- Reliable & durable
- Latex free
- Complete with oxygen tubing, reservoir bag & suitably sized anatomical air cushion mask
- Colour coded sizing guide on label
- Tamper proof, inner sealed polybag


DESCRIPTION	PRODUCT CODE	NPC CODE
1 Adult BVM Single Use 1500ml (Size 5 Mask)	TDM-RES-1500	FDE1629
2 Adult BVM Single Use 1500ml (Size 5 Mask w/Pop Off)	TDM-RES-1500W	FDE1625
Adult BVM Single Use 1000ml (Size 5 Mask)	TDM-RES-1000	FDE1628
3 Child BVM Single Use 550ml (Size 1 & 2 Masks w/Pop Off)	TDM-RES-0550	FDE1626
4 Infant BVM Single Use 280ml (Size 0 Mask w/Pop Off)	TDM-RES-0280	FDE1627



Adult 1500ml (Blue)



Adult 1000ml (Green)



Child 550ml (Pink)



Infant 280ml (Yellow)



SINGLE-USE AIR CUSHION ANAESTHETIC FACE MASK



- Available in neonate, infant, paediatric, small adult, medium adult, large adult & extra large adult
- The transparent mask allows for clear patient visibility
- Supplied without hook rings

DESCRIPTION	PRODUCT CODE	NPC CODE
Face Mask Size 0	TDM-MP-1600	FDD85370
Face Mask Size 1	TDM-MP-1601	FDD85371
Face Mask Size 2	TDM-MP-1602	FDD85374
Face Mask Size 3	TDM-MP-1603	FDD85375
Face Mask Size 4	TDM-MP-1604	FDD85372
Face Mask Size 5	TDM-MP-1605	FDD85373
Face Mask Size 6	TDM-MP-1606	FDD85377



SINGLE-USE AIR CUSHION ANAESTHETIC FACE MASK & CHECK VALVE



- Adjustable check valve ensuring minimum leakage & maximum patient comfort
- Available in neonate, infant, paediatric, small adult, medium adult, large adult & extra large adult
- The transparent mask allows for clear patient visibility
- Supplied without hook rings

DESCRIPTION	PRODUCT CODE	NPC CODE
Face Mask With Check Valve Size 0	TDM-MP-1610	FDD85378
Face Mask With Check Valve Size 1	TDM-MP-1611	FDD85376
Face Mask With Check Valve Size 2	TDM-MP-1612	FDD85382
Face Mask With Check Valve Size 3	TDM-MP-1613	FDD85383
Face Mask With Check Valve Size 4	TDM-MP-1614	FDD85380
Face Mask With Check Valve Size 5	TDM-MP-1615	FDD85381
Face Mask With Check Valve Size 6	TDM-MP-1616	FDD85379



Size 5

Size 0

REUSABLE SILICONE FACE MASKS

- Transparent to allow for visual checks
- Sizes range from premature to large adult
- 00, 0, 1 & 2 circular
- 3, 4, 5 anatomical
- Material: Silicone

DESCRIPTION	PRODUCT CODE	NPC CODE
Premature Mask Size 00	TDM-MP-1506	FDD5161
Neonatal Mask, Size 0	TDM-MP-1500	FDD668
Infant Mask Size 1	TDM-MP-1501	FDD669
Child Mask, Size 2	TDM-MP-1502	FDD670
Large Child/ Small Adult Mask, Size 3	TDM-MP-1513	FDD671
Medium Adult Mask, Size 4	TDM-MP-1514	FDD1055
Large Adult Mask, Size 5	TDM-MP-1515	FAG7008



TIMESCO OPTIMA GEUDAL OROPHARYNGEAL AIRWAY

STERILE | 2

- Latex Free & DEHP Free
- Individually wrapped
- Rounded atraumatic edges
- Bite blocks
- Available from size 000 (40mm) to size 6 (120mm)

COLOUR	DESCRIPTION	PRODUCT CODE	NPC CODE
PINK	Geudal Airway Size 000 (40mm)	TDM-D-4000	FEM85001
BLUE	Geudal Airway Size 00 (50mm)	TDM-D-400	FEM85000
BLACK	Geudal Airway Size 0 (60mm)	TDM-D-40	FEM85002
WHITE	Geudal Airway Size 1 (70mm)	TDM-D-4001	FEM85005
GREEN	Geudal Airway Size 2 (80mm)	TDM-D-4002	FEM85006
YELLOW	Geudal Airway Size 3 (90mm)	TDM-D-4003	FEM85003
RED	Geudal Airway Size 4 (100mm)	TDM-D-4004	FEM85004
ORANGE	Geudal Airway Size 5 (110mm)	TDM-D-4003	FEM85007
PURPLE	Geudal Airway Size 6 (120mm)	TDM-D-4004	FEM85008



NASOPHARYNGEAL AIRWAYS

STERILE | 2

- Extra large flange design, which means the device sits comfortably for the patient & there is no need for a safety pin
- Smooth surface permits easy insertion & withdrawal
- Single-use
- Latex-free
- Supplied in boxes of 10


DESCRIPTION	PRODUCT CODE	NPC CODE
Nasopharyngeal Airway Size 6	TDM-D-316	N/A
Nasopharyngeal Airway Size 7	TDM-D-317	N/A
Nasopharyngeal Airway Size 8	TDM-D-318	N/A
Nasopharyngeal Airway Size 9	TDM-D-319	N/A



BOUGIE INTRODUCERS (BOX OF 10)

STERILE | 2


DESCRIPTION	PRODUCT CODE	NPC CODE
Bougie Straight 5Ch 47cm Neonate	TDM-D-1130	FSM742
Bougie Straight 10Ch 60cm Child	TDM-D-1131	FSM743
Bougie Angled 15Ch 80cm Adult	TDM-D-1138	FSM744
Flexible Tip Boujie 65cm	TDM-D-1150	N/A



ECO MAGILL FORCEPS (BOX OF 10)

STERILE | 2

DESCRIPTION	PRODUCT CODE	NPC CODE
Eco Magill Forceps Infant 16cm	D.04.200.15	N/A
Eco Magill Forceps Child 20cm	D.04.200.25	N/A
Eco Magill Forceps Adult 24.5cm	D.04.200.35	N/A



ENDOTRACHEAL TUBES

STERILE 2

- Available in sizes 3.0 – 10.0 (cuffed)
- Available in sizes 2.0 – 5.0 (uncuffed)
- Includes Murphy Eye
- High quality inflation control valve (cuffed)
- Clear 15mm connector
- Radio-opaque line for identification on X-Ray
- Depth marker lines
- Smooth soft tip reduces the potential of tracheal trauma to a patient
- Supplied in boxes of 20
- Latex Free

DESCRIPTION	PRODUCT CODE	NPC CODE
Endotracheal Tube 2.0 (x20)	TDM-D-5220	N/A
Endotracheal Tube 2.5 (x20)	TDM-D-5225	N/A
Endotracheal Tube 3.0 (x20)	TDM-D-5230	N/A
Endotracheal Tube 3.5 (x20)	TDM-D-5235	N/A
Endotracheal Tube 4.0 (x20)	TDM-D-5240	N/A
Endotracheal Tube 4.5 (x20)	TDM-D-5245	N/A
Endotracheal Tube 5.0 (x20)	TDM-D-5250	N/A
Endotracheal Tube Cuffed 3.0 (x20)	TDM-D-5130	FSM85000
Endotracheal Tube Cuffed 3.5 (x20)	TDM-D-5135	FSM85002
Endotracheal Tube Cuffed 4.0 (x20)	TDM-D-5140	FSM85008
Endotracheal Tube Cuffed 4.5 (x20)	TDM-D-5145	FSM85010

DESCRIPTION	PRODUCT CODE	NPC CODE
Endotracheal Tube Cuffed 5.0 (x20)	TDM-D-5150	FSM85006
Endotracheal Tube Cuffed 5.5 (x20)	TDM-D-5155	FSM85007
Endotracheal Tube Cuffed 6.0 (x20)	TDM-D-5160	FSM85013
Endotracheal Tube Cuffed 6.5 (x20)	TDM-D-5165	FSM85014
Endotracheal Tube Cuffed 7.0 (x20)	TDM-D-5170	FSM85012
Endotracheal Tube Cuffed 7.5 (x20)	TDM-D-5175	FSM85004
Endotracheal Tube Cuffed 8.0 (x20)	TDM-D-5180	FSM85005
Endotracheal Tube Cuffed 8.5 (x20)	TDM-D-5185	FSM85001
Endotracheal Tube Cuffed 9.0 (x20)	TDM-D-5190	FSM85003
Endotracheal Tube Cuffed 9.5 (x20)	TDM-D-5195	FSM85009
Endotracheal Tube Cuffed 10.0 (x20)	TDM-D-5199	FSM85011



DISPOSABLE SILICONE LARYNGEAL MASK

STERILE 2

- Low pressure seal around laryngeal inlet, permitting gentle positive pressure ventilation
- Sterile, Latex Free & DEHP Free
- Medical grade silicone

DESCRIPTION	PRODUCT CODE	NPC CODE
Laryngeal Mask Silicone SU (Size 1)	TDM-LMO-110	FXU35042
Laryngeal Mask Silicone SU (Size 1.5)	TDM-LMO-115	FXU35037
Laryngeal Mask Silicone SU (Size 2)	TDM-LMO-120	FXU35041
Laryngeal Mask Silicone SU (Size 2.5)	TDM-LMO-125	FXU35045
Laryngeal Mask Silicone SU (Size 3)	TDM-LMO-130	FXU35040
Laryngeal Mask Silicone SU (Size 4)	TDM-LMO-140	FXU35038
Laryngeal Mask Silicone SU (Size 5)	TDM-LMO-150	FXU35039



CALLISTO PRELOADED HANDLE



- Powerful, bright LED bulb for excellent illumination
- Low volume plastic construction reducing plastic waste & cost of disposal
- Single use
- Preloaded with alkaline button cell batteries quick & easy to remove for safe battery disposal
- Green ISO 7376 fitting compatible with green fitting laryngoscope blades
- Metal hook bar offering exceptional strength for a single use handle
- Supplied in boxes of 25

DESCRIPTION	PRODUCT CODE	NPC CODE
Callisto Preloaded Handle, Compact	DS.3940.460.10	FSM3228
Callisto Preloaded Handle, V2	DS.3940.450.10	FSM85017



XENON LARYNGOSCOPE HANDLE

- A polished metal collar reflects and directs light for maximum illumination
- High performance replaceable 2.5V Xenon bulb
- The Xenon handle works with dry cell or rechargeable batteries. Secure screw cap for easy access to batteries
- Knurled handle design provides excellent grip
- Colour coded band to clearly identify it is ISO 7376 fitting compatible with green fitting blades
- Clearly marked with the Xenon logo, product code, CE mark and LOT number.
- Fully autoclavable (excluding bulb)

DESCRIPTION	PRODUCT CODE	NPC CODE
Xenon Laryngoscope Handle, Paediatric	3000.350.05	FSM426
Xenon Laryngoscope Handle, Adult	3000.350.10	FSM065



XLED LARYNGOSCOPE HANDLE

- A polished metal collar reflects and directs light for maximum illumination
- High performance replaceable 3v LED bulb
- The XLED handle works with dry cell or rechargeable batteries. Secure screw cap for easy access to batteries
- Knurled handle design provides excellent grip.
- Colour coded band to clearly identify it is ISO 7376 fitting compatible with green fitting blades
- Clearly marked with the XLED logo, product code, CE mark and LOT number.
- Fully autoclavable (excluding bulb)

DESCRIPTION	PRODUCT CODE	NPC CODE
Xenon Laryngoscope Handle, Paediatric	3000.500.10	FSM717
Xenon Laryngoscope Handle, Adult	3000.500.05	FSM718

CALLISTO BLADES



- CE marked
- Single-Use
- Latex Free
- ISO 7376 compatible
- Supplied with multilingual instructions
- Supplied in boxes of 10



DESCRIPTION	PRODUCT CODE	NPC CODE
Callisto Laryn Blade Macintosh 1	DS.3940.150.10	FSM509
Callisto Laryn Blade Macintosh 2	DS.3940.150.15	FSM076
Callisto Laryn Blade Macintosh 3	DS.3940.150.20	FSM077
Callisto Laryn Blade Macintosh 4	DS.3940.150.25	FSM078
Callisto Laryn Blade Miller 000	DS.3940.185.02	FSM4247
Callisto Laryn Blade Miller 00	DS.3940.185.03	FSM763
Callisto Laryn Blade Miller 0	DS.3940.185.05	FSM103
Callisto Laryn Blade Miller 1	DS.3940.185.10	FSM104
Callisto Laryn Blade Miller 2	DS.3940.185.15	FSM106
Callisto Laryn Blade Miller 3	DS.3940.185.20	FSM107

CALLISTO BLADE & PRELOADED COMPACT HANDLE SET



- Powerful, bright LED bulb for excellent illumination
- Handle measurements: 139mm (H) x 33mm (Ø)
- Preloaded with alkaline button cell batteries
Quick and easy to remove for safe battery disposal
- Supplied in boxes of 10
- Metal hook bar offering exceptional strength for a single use handle
- Green ISO 7376 fitting compatible with green fitting laryngoscope blades
- Single use



DESCRIPTION	PRODUCT CODE	NPC CODE
Callisto Macintosh Blade Size 2 & Callisto Preloaded Compact Handle Set	DS.3950.952.15	FSM3229
Callisto Macintosh Blade Size 3 & Callisto Preloaded Compact Handle Set	DS.3950.952.20	FSM3230
Callisto Macintosh Blade Size 4 & Callisto Preloaded Compact Handle Set	DS.3950.952.25	FSM3231
Callisto Miller Blade Size 00 & Callisto Preloaded Compact Handle Set	DS.3950.952.28	FSM3814
Callisto Miller Blade Size 0 & Callisto Preloaded Compact Handle Set	DS.3950.952.30	FSM3232
Callisto Miller Blade Size 1 & Callisto Preloaded Compact Handle Set	DS.3950.952.35	FSM3233
Callisto Miller Blade Size 2 & Callisto Preloaded Compact Handle Set	DS.3950.952.40	FSM3815
Callisto Miller Blade Size 3 & Callisto Preloaded Compact Handle Set	DS.3950.952.45	FSM3816



OPTIMA VIEW VIDEO LARYNGOSCOPE

FOR USE WITH OPTIMA VIEW SINGLE USE BLADES

- Rechargeable battery lasting for over 2 hours
- HDMI connectivity to external monitor or TV screen
- The screen rotates 180° from left to right & tilts 110° up & down. It also has a view angle of 160° making it convenient for teaching with multiple viewers
- Blades are available in sizes Mac 1 – Mac 5 including channelled Mac 3 – Mac 5. Shaped similarly to traditional Macintosh blades
- The 2.0 megapixel, high resolution camera features a 50°- 60° lens angle
- Resolution ratio: 800 x 480RGB
- Built-in 7.5GB storage, which can store approx. 10 hours of video or 70k+ pictures
- One-click button suitably placed on the ergonomic handle for photography & video recording
- The sterile, single-use blade simply slides & clips onto the Optima View system

DESCRIPTION	PRODUCT CODE	NPC CODE
Optima View Video SU Lary Handle	4000.310.50	FSM85626
Optima View Video Lary Carry Case	4000.900.13	N/A




OPTIMA VIEW VIDEO LARYNGOSCOPE

FOR USE WITH OPTIMA VIEW REUSABLE BLADES

- Rechargeable battery lasting for over 2 hours
- HDMI connectivity to external monitor or TV screen
- The screen rotates 180° from left to right & tilts 110° up & down. It also has a view angle of 160° making it convenient for teaching with multiple viewers
- Blades are available in sizes Mac 1 – Mac 5 including channelled Mac 3 – Mac 5. Shaped similarly to traditional Macintosh blades
- The 2.0 megapixel, high resolution camera features a 50°- 60° lens angle
- Resolution ratio: 800 x 480RGB
- Built-in 7.5GB storage, which can store approx. 10 hours of video or 70k+ pictures
- One-click button suitably placed on the ergonomic handle for photography & video recording
- Traditional Hook bar system offering exceptional strength.

DESCRIPTION	PRODUCT CODE	NPC CODE
Optima View Video RU Lary Handle	4000.300.50	FSM85639
Optima View Video Lary Carry Case	4000.900.13	N/A




OPTIMA VIEW VIDEO LARYNGOSCOPE BLADE SINGLE-USE

STERILE

- Compatible with Optima View Video Laryngoscope handle
- Supplied with instructions
- Polycarbonate
- CE marked

	DESCRIPTION	PRODUCT CODE	NPC CODE
1	Optima View Video Lary Single-use Mac 1 Blade, Pack of 10	DS.4500.150.10X	FSM85625
2	Optima View Video Lary Single-use Mac 2 Blade, Pack of 10	DS.4500.150.15X	FSM85635
3	Optima View Video Lary Single-use Mac 3 Blade, Pack of 10	DS.4500.150.20X	FSM85638
4	Optima View Video Lary Single-use Mac 4 Blade, Pack of 10	DS.4500.150.25X	FSM85633
5	Optima View Video Lary Single-use Mac 5 Blade, Pack of 10	DS.4500.150.30X	N/A



OPTIMA VIEW VIDEO LARYNGOSCOPE BLADE REUSABLE

- Compatible with Optima View Video Laryngoscope handle
- Supplied with instructions
- Stainless Steel
- CE marked

	DESCRIPTION	PRODUCT CODE	NPC CODE
6	Optima View Video Lary Reusable Mac 1 Blade, Single	4500.150.10	FSM85630
7	Optima View Video Lary Reusable Mac 2 Blade, Single	4500.150.15	FSM85623
8	Optima View Video Lary Reusable Mac 3 Blade, Single	4500.150.20	FSM85624
9	Optima View Video Lary Reusable Mac 4 Blade, Single	4500.150.25	FSM85622
10	Optima View Video Lary Reusable Mac 5 Blade, Single	4500.150.30	N/A
11	Optima View Video Lary Reusable Miller 00 Blade, Single	4500.185.03	FSM85634
12	Optima View Video Lary Reusable Miller 0 Blade, Single	4500.185.05	FSM85643
13	Optima View Video Lary Reusable Miller 1 Blade, Single	4500.185.10	FSM85644

CHEST DRAINAGE UNIT



★ NEW
PRODUCT

‘WAVE’ 3-CHAMBER WATER SEAL DRAIN

STERILE | 2

The Vide® Wave Water Seal Chest Drainage Unit is a 2000ml, three-chamber chest drain featuring:
Three-column collection chamber
A water seal/leak monitor with a suction control chamber.
Providing precise and exact suction control up to a maximum vacuum level of 25cmH2O.

The Wave chest drain integrates a variety of safety features designed to protect patient from sudden changes in the intrapleural environment.

DESCRIPTION	PRODUCT CODE	NPC CODE
‘Wave’ 3-Chamber Water Seal Chest Drain - 2ltr (box of 5)	D50.100	FDR85009



★ NEW
PRODUCT


‘WAVE’ PAEDIATRIC WATER SEAL DRAIN

STERILE | 2

The Vide® Paediatric Wave Water Seal Chest Drainage Unit is a 200ml, three-chamber chest drain featuring:
Three-column collection chamber
A water seal/leak monitor with a suction control chamber.
Providing precise and exact suction control up to a maximum vacuum level of 25cmH2O.

The Wave chest drain integrates a variety of safety features designed to protect patient from sudden changes in the intrapleural environment.

DESCRIPTION	PRODUCT CODE	NPC CODE
‘Wave’ Paediatric Water Seal Chest Drain - 200ml (box of 5)	D50.115	FDR85008



★ NEW
PRODUCT

'SHORE' DRY SUCTION / WATER SEAL DRAIN

STERILE

2

The Vide® Shore chest drain is a 2000ml, three-chamber chest drain unit that offers the reliability of a traditional water seal drain and the convenience of a water-less suction control operation.

The Shore chest drain completely eliminates the need for constantly maintaining water level in the suction control chamber (and thus the need to scramble for sterile water in the midst of a night!). The elimination of water also means the vacuum level can be maintained with less fluctuations than a traditional water-based suction control.

DESCRIPTION	PRODUCT CODE	NPC CODE
'Shore' Dry Suction / Water Seal Chest Drain (box of 5)	D50.110	N/A



★ NEW
PRODUCT

'SAND' DRY SEAL DRAIN

STERILE

2

The Vide® Sand Chest Drain is a 2000ml three-chamber drain that offers all of the benefits and features of our traditional drains but requires minimal manual intervention, waterless technology and enhanced safety reduces the staff burden.

The automatic high negative pressure relief valve limits system pressure.

Knock-over nozzles allow fluids levels in the air leak monitor chamber to be recovered instantly.

Positive relief valve is mounted to ensure pressure build-up can be relieved even during external suction disruption.

DESCRIPTION	PRODUCT CODE	NPC CODE
'Sand' Dry Seal Chest Drain (box of 5)	D50.105	FDR85010

GO GREEN.

We accept that our daily activities may have both a positive & negative impact on the environment. We are committed to managing our operations, products & services in such a way to preserve, protect & enhance the environment. Our goal is to minimise this effect & to reduce our carbon footprint wherever possible.

We aim to achieve this by:

- Ensuring that every aspect of our activities is conducted in accordance with environmental policies.
- Minimising the consumption of natural resources & energy, whilst consuming material goods in moderation.
- Reducing the creation of waste by improving practices & recycling materials whenever practical.
- Continually seeking to improve environmental practices & performance & develop more environmentally friendly products.
- Informing & educating all Timesco staff about environmental issues in the context of the business through appropriate training & information sharing.

Timesco Healthcare Ltd., 3 Carnival Park, Carnival Close,
Basildon, Essex, SS14 3WN, United Kingdom
T: +44 (0)1268 297 700 E: info@timesco.com W: timesco.com



TIMESCO