



ANAESTHESIA  
TIMESCO



**AIRWAY  
MANAGEMENT**  
TIMESCO™

<b>PRODUCT RANGE</b>	<b>PAGE</b>
Nebulisers / Respiratory	4
Resuscitation Aids	5
Face Masks	6
Airway Management	7 - 8
Magill Forceps	9



★ NEW PRODUCT

## NEBULISER KIT (MOUTHPIECE)

- Supplied with; mouthpiece, O2 tubing & nebuliser medication chamber (6ml)
- Light weight medical grade PVC
- Latex & Phthalates free
- Box of 50

DESCRIPTION	PRODUCT CODE	NPC CODE
Nebuliser Kit with Mouthpiece	D10.017	FAG85148



★ NEW PRODUCT

## ADULT NON-REBREATH MASK

- Supplied with; mask, O2 tubing & reservoir bag
- With nose clip & elastic head strap
- Light weight medical grade PVC
- Latex & Phthalates free
- High concentration 99% to 100% delivery
- Box of 50

DESCRIPTION	PRODUCT CODE	NPC CODE
Adult Non-Rebreath Mask	D10.305	FAG7028



★ NEW PRODUCT

## PAEDIATRIC NON-REBREATH MASK

- Supplied with; mask, O2 tubing & reservoir bag
- With nose clip & elastic head strap
- Light weight medical grade PVC
- Latex & Phthalates free
- High concentration 99% to 100% delivery
- Box of 50

DESCRIPTION	PRODUCT CODE	NPC CODE
Paediatric Non-Rebreath Mask	D10.306	FAG7025



## CPR RESCU-MASK



- Single-use
- Replacement filters available for training purposes

Supplied with:

- One-way bacterial & viral filter
- Built-in O2 port
- Elasticated head strap
- Hard plastic carry case

DESCRIPTION	PRODUCT CODE	NPC CODE
CPR Rescu-Mask	TDM-PM-250	FDD2552
One-way valve / filter (x5)	TDM-PM-110	FAG7013



## SINGLE-USE RESUSCITATORS (BVM)



- Textured grip (non-slip)
- 360° swivel connector for greater versatility
- Glued components (with exception of face mask)
- Adjustable air cushion face mask for maximum patient comfort
- Transparent bag for complete internal visibility
- Excellent recoil recovery following bag compression
- Available in a variety of sizes
- Reliable & durable
- Latex free
- Complete with oxygen tubing, reservoir bag & suitably sized anatomical air cushion mask
- Colour coded sizing guide on label
- Tamper proof, inner sealed polybag


DESCRIPTION	PRODUCT CODE	NPC CODE
1 Adult BVM Single Use 1500ml (Size 5 Mask)	TDM-RES-1500	FDE1629
2 Adult BVM Single Use 1500ml (Size 5 Mask w/Pop Off)	TDM-RES-1500W	FDE1625
Adult BVM Single Use 1000ml (Size 5 Mask)	TDM-RES-1000	FDE1628
3 Child BVM Single Use 550ml (Size 1 & 2 Masks w/Pop Off)	TDM-RES-0550	FDE1626
4 Infant BVM Single Use 280ml (Size 0 Mask w/Pop Off)	TDM-RES-0280	FDE1627




Adult 1500ml (Blue)



Adult 1000ml (Green)



Child 550ml (Pink)



Infant 280ml (Yellow)



## SINGLE-USE AIR CUSHION ANAESTHETIC FACE MASK



- Available in neonate, infant, paediatric, small adult, medium adult, large adult & extra large adult
- The transparent mask allows for clear patient visibility
- Supplied without hook rings

DESCRIPTION	PRODUCT CODE	NPC CODE
Face Mask Size 0	TDM-MP-1600	<a href="#">FAG7016</a>
Face Mask Size 1	TDM-MP-1601	<a href="#">FAG7019</a>
Face Mask Size 2	TDM-MP-1602	<a href="#">FAG7023</a>
Face Mask Size 3	TDM-MP-1603	<a href="#">FAG7009</a>
Face Mask Size 4	TDM-MP-1604	<a href="#">FAG7034</a>
Face Mask Size 5	TDM-MP-1605	<a href="#">FAG7017</a>
Face Mask Size 6	TDM-MP-1606	<a href="#">FAG7029</a>



## SINGLE-USE AIR CUSHION ANAESTHETIC FACE MASK & CHECK VALVE



- Adjustable check valve ensuring minimum leakage & maximum patient comfort
- Available in neonate, infant, paediatric, small adult, medium adult, large adult & extra large adult
- The transparent mask allows for clear patient visibility
- Supplied without hook rings

DESCRIPTION	PRODUCT CODE	NPC CODE
Face Mask With Check Valve Size 0	TDM-MP-1610	<a href="#">FAG7033</a>
Face Mask With Check Valve Size 1	TDM-MP-1611	<a href="#">FAG7032</a>
Face Mask With Check Valve Size 2	TDM-MP-1612	<a href="#">FAG7027</a>
Face Mask With Check Valve Size 3	TDM-MP-1613	<a href="#">FAG7014</a>
Face Mask With Check Valve Size 4	TDM-MP-1614	<a href="#">FAG7030</a>
Face Mask With Check Valve Size 5	TDM-MP-1615	<a href="#">FAG7020</a>
Face Mask With Check Valve Size 6	TDM-MP-1616	<a href="#">FAG7012</a>



Size 5

Size 0

## REUSABLE SILICONE FACE MASKS

- Transparent to allow for visual checks
- Sizes range from premature to large adult
- 00, 0, 1 & 2 circular
- 3, 4, 5 anatomical
- Material: Silicone

DESCRIPTION	PRODUCT CODE	NPC CODE
Premature Mask Size 00	TDM-MP-1506	<a href="#">FDD5161</a>
Neonatal Mask, Size 0	TDM-MP-1500	<a href="#">FDD668</a>
Infant Mask Size 1	TDM-MP-1501	<a href="#">FDD669</a>
Child Mask, Size 2	TDM-MP-1502	<a href="#">FDD670</a>
Large Child/ Small Adult Mask, Size 3	TDM-MP-1513	<a href="#">FDD671</a>
Medium Adult Mask, Size 4	TDM-MP-1514	<a href="#">FDD1055</a>
Large Adult Mask, Size 5	TDM-MP-1515	<a href="#">FAG7008</a>




## GUEDEL AIRWAYS



- Latex Free
- Individually wrapped
- Rounded atraumatic edges
- Bite blocks
- Available from size 000 to size 4

DESCRIPTION	PRODUCT CODE	NPC CODE
Guedal Airway Size 000	TDM-D-3190	N/A
Guedal Airway Size 00	TDM-D-3195	N/A
Guedal Airway Size 0	TDM-D-3200	N/A
Guedal Airway Size 1	TDM-D-3201	N/A
Guedal Airway Size 2	TDM-D-3202	N/A
Guedal Airway Size 3	TDM-D-3203	N/A
Guedal Airway Size 4	TDM-D-3204	N/A



## NASOPHARYNGEAL AIRWAYS



- Extra large flange design, which means the device sits comfortably for the patient & there is no need for a safety pin
- The one piece design is made from a transparent material, which enables blockages to be easily & quickly identified
- Single-use
- Latex-free
- Supplied in boxes of 10

DESCRIPTION	PRODUCT CODE	NPC CODE
Nasopharyngeal Airway Size 6	TDM-D-106	FSM4261
Nasopharyngeal Airway Size 7	TDM-D-107	FSM4259
Nasopharyngeal Airway Size 8	TDM-D-108	FSM4260
Nasopharyngeal Airway Size 9	TDM-D-109	FSM4340




## BOUGIE INTRODUCERS



- Bougie introducers will assist in directional control of the endotracheal tube, especially when the laryngeal opening is not fully visible
- Latex-free
- Sterile and individually wrapped
- Available in neonate, child and adult sizes
- Single-use
- Supplied in boxes of 10

DESCRIPTION	PRODUCT CODE	NPC CODE
Bougie Straight 5Ch 470mm Neonate (Box of 10)	TDM-D-1130	FSM742
Bougie Straight 10Ch 600mm Child (Box of 10)	TDM-D-1131	FSM743
Bougie Angled 15Ch 800mm Adult (Box of 10)	TDM-D-1138	FSM744



## ENDOTRACHEAL TUBES



- Available in sizes 3.0 – 10.0 (cuffed)
- Includes Murphy Eye
- High quality inflation control valve
- Clear 15mm connector
- Radio-opaque line for identification on X-Ray
- Single Use
- Depth marker lines
- High volume, low pressure cuff
- Smooth soft tip reduces the potential of tracheal trauma to a patient
- Supplied in boxes of 20

DESCRIPTION	PRODUCT CODE	NPC CODE
Endotracheal Tube Cuffed 3.0 (x20)	TDM-D-5130	<a href="#">FSM85000</a>
Endotracheal Tube Cuffed 3.5 (x20)	TDM-D-5135	<a href="#">FSM85002</a>
Endotracheal Tube Cuffed 4.0 (x20)	TDM-D-5140	<a href="#">FSM85008</a>
Endotracheal Tube Cuffed 4.5 (x20)	TDM-D-5145	<a href="#">FSM85010</a>
Endotracheal Tube Cuffed 5.0 (x20)	TDM-D-5150	<a href="#">FSM85006</a>
Endotracheal Tube Cuffed 5.5 (x20)	TDM-D-5155	<a href="#">FSM85007</a>
Endotracheal Tube Cuffed 6.0 (x20)	TDM-D-5160	<a href="#">FSM85013</a>
Endotracheal Tube Cuffed 6.5 (x20)	TDM-D-5165	<a href="#">FSM85014</a>
Endotracheal Tube Cuffed 7.0 (x20)	TDM-D-5170	<a href="#">FSM85012</a>
Endotracheal Tube Cuffed 7.5 (x20)	TDM-D-5175	<a href="#">FSM85004</a>
Endotracheal Tube Cuffed 8.0 (x20)	TDM-D-5180	<a href="#">FSM85005</a>
Endotracheal Tube Cuffed 8.5 (x20)	TDM-D-5185	<a href="#">FSM85001</a>
Endotracheal Tube Cuffed 9.0 (x20)	TDM-D-5190	<a href="#">FSM85003</a>
Endotracheal Tube Cuffed 9.5 (x20)	TDM-D-5195	<a href="#">FSM85009</a>
Endotracheal Tube Cuffed 10.0 (x20)	TDM-D-5199	<a href="#">FSM85011</a>


## DISPOSABLE SILICONE LARYNGEAL MASK



- Low pressure seal around laryngeal inlet, permitting gentle positive pressure ventilation
- Sterile device
- Medical grade silicone

DESCRIPTION	PRODUCT CODE	NPC CODE
Laryngeal Mask Silicone SU (Size 1)	TDM-LMO-110	<a href="#">FSM3160</a>
Laryngeal Mask Silicone SU (Size 1.5)	TDM-LMO-115	<a href="#">FSM3161</a>
Laryngeal Mask Silicone SU (Size 2)	TDM-LMO-120	<a href="#">FSM3162</a>
Laryngeal Mask Silicone SU (Size 2.5)	TDM-LMO-125	<a href="#">FSM3163</a>
Laryngeal Mask Silicone SU (Size 3)	TDM-LMO-130	<a href="#">FSM3164</a>
Laryngeal Mask Silicone SU (Size 4)	TDM-LMO-140	<a href="#">FSM3165</a>
Laryngeal Mask Silicone SU (Size 5)	TDM-LMO-150	<a href="#">FSM3166</a>


## SINGLE-USE ECO MAGILL FORCEPS



- CE Marked
- Latex Free
- Place/remove throat packs and clear foreign bodies from the upper airway.
- Stainless steel with the look and feel of reusable counterparts.
- Angled to maneuver the ET tube in nasal intubations.
- Available in 3 sizes to meet the specific requirements of the patients.
- Single use to reduce the risk of cross infection.
- Eye design for easy ET tube maneuverability and serrated tip for excellent grip.
- Box of 10

DESCRIPTION	PRODUCT CODE	NPC CODE
Eco Magill Forceps Infant 16cm	D.04.200.15	N/A
Eco Magill Forceps Child 20cm	D.04.200.25	N/A
Eco Magill Forceps Adult 24.5cm	D.04.200.35	N/A

## REUSABLE MAGILL FORCEPS



- CE Marked
- Latex Free
- Place/remove throat packs and clear foreign bodies from the upper airway.
- Stainless steel with the look and feel of reusable counterparts.
- Angled to maneuver the ET tube in nasal intubations.
- Available in 3 sizes to meet the specific requirements of the patients.
- Eye design for easy ET tube maneuverability and serrated tip for excellent grip.

DESCRIPTION	PRODUCT CODE	NPC CODE
Reusable Magill Forceps Infant 15.5cm	04.200.10	FFQ4592
Reusable Magill Forceps Child 18.5cm	04.200.20	FFQ4593
Reusable Magill Forceps Adult 23cm	04.200.30	FFQ4594



Timesco Healthcare Ltd., 3 Carnival Park, Carnival Close, Basildon, Essex, SS14 3WN, United Kingdom  
T: +44 (0)1268 297 700 E: [info@timesco.com](mailto:info@timesco.com) W: [timesco.com](http://timesco.com)

B.GM.015/3 AUG25 E & O.E. Timesco Healthcare Ltd.